

# Laboratory Chairs and Stools



# CHAIRS AND STOOLS

Our products are manufactured by hand to ensure quality that will last. Our chairs' durability is suitable for even the harshest of environments, and each chair comes with a manufacturer's warranty.

## Features include:

- **Ergonomics:** Easyjust controls put everything within intuitive reach. Adjustable footrests help circulation and ease fatigue. Optional armrests aid in preventing cumulative trauma disorders, while easy height adjustments put anyone on proper eye level.
- **Materials:** Our products are made for long life, with a choice in bases of heavy-duty chromed aluminum or nylon reinforced fiberglass. Seat and back materials include customizable fabric or vinyl, super-tough polyurethane, or seven-ply hardwood.
- **Colors:** Wide array of fabric coverings that will complement any décor, and offers 100's of colors of vinyl to select from, with seven standard colors for fabric and vinyl.
- **Modularity:** High quality interlocking and interchangeable product components provides for a competitive edge by creating a customizable configuration array that answers any type of workplace seating demands and allows for rapid order turnaround.

## Explore different Categories

- **Laboratory Chairs and Stools**
- **Office chairs and Stools**
- **Industrial Chair and Stools**
- **Medical chairs and Stools**
- **Clean room Chairs and stools**



# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS- VDHCH-RG-T0-A0-RC-8613
<b>DESCRIPTION</b>	Vinyl Chair - Desk Height with Casters in Orange Kist Trailblazer™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18-23 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 2.5in
<b>BACK</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>AVAILABLE COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

<b>CATALOG NO.</b>	MLS-VMBCH-RG-T0-A0-CF-RC-8613
<b>DESCRIPTION</b>	Vinyl Chair - Medium Bench Height with Chrome Foot Ring & Casters in Orange Kist Trailblazer™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21-28.25 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 2.5in
<b>BACK</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>AVAILABLE COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



<b>CATALOG NO.</b>	MLS-VHBCH-RG-T0-A0-CF-RC-8613
<b>DESCRIPTION</b>	Vinyl Chair - High Bench Height with Chrome Foot Ring & Casters in Orange Kist Trailblazer™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 25-34.5 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 2.5in
<b>BACK</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>AVAILABLE COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS- VDHCH-CR-T0-A0-CC-8801
<b>DESCRIPTION</b>	Vinyl Chair Chrome - Desk Height with Casters in Marine Blue Supernova™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19-24 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 2.5in
<b>BACK</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>AVAILABLE COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

<b>CATALOG NO.</b>	MLS- VMBCH-CR-T0-A0-CF-CC-8801
<b>DESCRIPTION</b>	Vinyl Chair Chrome - Medium Bench Height with Chrome Foot Ring & Casters in Marine Blue Supernova™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21.5-29 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 2.5in
<b>BACK</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>AVAILABLE COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



<b>CATALOG NO.</b>	MLS- VDHCH-CR-T0-A0-CC-8801
<b>DESCRIPTION</b>	Vinyl Chair Chrome - High Bench Height with Chrome Foot Ring & Casters in Marine Blue Supernova™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 26.25-36 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 2.5in
<b>BACK</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>AVAILABLE COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-VDHCH-RT-T0-A0-RG-8584
<b>DESCRIPTION</b>	Vinyl Chair - Desk Height with Round Tube Base & Stationary Glides in Taupe Trailblazer™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.5-25.75 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 2.5in
<b>BACK</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>AVAILABLE COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

<b>CATALOG NO.</b>	MLS-VMBCH-RT-T0-A0-RG-8584
<b>DESCRIPTION</b>	Vinyl Chair - Medium Bench Height with Round Tube Base & Stationary Glides in Taupe Trailblazer™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 23.25-31 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 2.5in
<b>BACK</b>	Width: 16.5in, Height: 11.5in, Thickness: 2.5in
<b>AVAILABLE COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



<b>CATALOG NO.</b>	MLS-WDHCH-RG-RC
<b>DESCRIPTION</b>	Wood Chair - Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 15-20.5 in
<b>SEAT</b>	Width: 17.25in, Depth: 15.5 in, Thickness: 0.5in
<b>BACK</b>	Width: 15.5in, Height: 10in, Thickness: 0.5in
<b>COLOR</b>	Blonde Wood



# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-WMBCH-RG-CF-RC
<b>DESCRIPTION</b>	Wood Chair - Medium Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18-25.5 in
<b>SEAT</b>	Width: 17.25in, Depth: 15.5 in, Thickness: 0.5in
<b>BACK</b>	Width: 15.5in, Height: 10in, Thickness: 0.5in
<b>COLOR</b>	Blonde Wood

<b>CATALOG NO.</b>	MLS-WHBCH-RG-CF-RC
<b>DESCRIPTION</b>	Wood Chair - High Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 22-32 in
<b>SEAT</b>	Width: 17.25in, Depth: 15.5 in, Thickness: 0.5in
<b>BACK</b>	Width: 15.5in, Height: 10in, Thickness: 0.5in
<b>COLOR</b>	Blonde Wood



<b>CATALOG NO.</b>	MLS-WDHCH-RT-RG
<b>DESCRIPTION</b>	Wood Chair - Desk Height with Round Tube Base & Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17.75-23 in
<b>SEAT</b>	Width: 17.25in, Depth: 15.5 in, Thickness: 0.5in
<b>BACK</b>	Width: 15.5in, Height: 10in, Thickness: 0.5in
<b>COLOR</b>	Blonde Wood

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS- WMBCH-RT-RG
<b>DESCRIPTION</b>	Wood Chair - Medium Bench Height with Round Tube Base & Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.5-28 in
<b>SEAT</b>	Width: 17.25in, Depth: 15.5 in, Thickness: 0.5in
<b>BACK</b>	Width: 15.5in, Height: 10in, Thickness: 0.5in
<b>COLOR</b>	Blonde Wood

<b>CATALOG NO.</b>	MLS-VDHCH-MB-RG-T0-A0-RC-8801
<b>DESCRIPTION</b>	Vinyl Chair - Desk Height with Medium Back & Casters in Marine Blue Supernova™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18-23 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 3in
<b>BACK</b>	Width: 16.5in, Height: 18in, Thickness: 2.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



<b>CATALOG NO.</b>	MLS-VMBCH-MB-RG-T0-A0-CF-RC-8801
<b>DESCRIPTION</b>	Vinyl Chair - Medium Bench Height with Medium Back, Chrome Foot Ring, & Casters in Marine Blue Supernova™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21-28.25 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 3in
<b>BACK</b>	Width: 16.5in, Height: 18in, Thickness: 2.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS- VDHCH-MB-CR-T0-A0-CC-8567
<b>DESCRIPTION</b>	Vinyl Chair Chrome – Desk Height with Medium Back & Casters in Dove Trailblazer™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19–24 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 3in
<b>BACK</b>	Width: 16.5in, Height: 18in, Thickness: 2.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

<b>CATALOG NO.</b>	MLS- VHBCH-MB-CR-T0-A0-CF-CC-8567
<b>DESCRIPTION</b>	Vinyl Chair Chrome – High Bench Height with Medium Back, Chrome Foot Ring, & Casters in Dove Trailblazer™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 26.25–36 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 3in
<b>BACK</b>	Width: 16.5in, Height: 18in, Thickness: 2.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



<b>CATALOG NO.</b>	MLS-VDHCH-MB-RT-T0-A0-RG-8567
<b>DESCRIPTION</b>	Vinyl Chair – Desk Height with Round Tube Base, Medium Back, & Stationary Glides in Dove Trailblazer™ Vinyl
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.5–25.75 in
<b>SEAT</b>	Width: 18.5in, Depth: 17in, Thickness: 3in
<b>BACK</b>	Width: 16.5in, Height: 18in, Thickness: 2.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-BPDHCH-RG-RC
<b>DESCRIPTION</b>	Core Polyurethane Chair – Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 15.5-20.75 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in

<b>CATALOG NO.</b>	MLS-BPMBCH-RG-CF-RC
<b>DESCRIPTION</b>	Core Polyurethane Chair – Medium Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18.25-26 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in



<b>CATALOG NO.</b>	MLS- BPHBCH-RG-CF-RC
<b>DESCRIPTION</b>	Core Polyurethane Chair – High Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 22.5-32.5 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in

<b>CATALOG NO.</b>	MLS-BPDHCH-CR-CC
<b>DESCRIPTION</b>	Core Polyurethane Chair Chrome – Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 16.75-22 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in



# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-BPMBCH-CR-CF-CC
<b>DESCRIPTION</b>	Core Polyurethane Chair Chrome – Medium Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19.5-27 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in

<b>CATALOG NO.</b>	MLS-BPHBCH-CR-CF-CC
<b>DESCRIPTION</b>	Core Polyurethane Chair Chrome – High Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 23.75-33.75 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in



<b>CATALOG NO.</b>	MLS-BPDHCH-RT-RG
<b>DESCRIPTION</b>	Core Polyurethane Chair – Desk Height with Round Tube Base and Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18-23.5 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in

<b>CATALOG NO.</b>	MLS- BPMBCH-RT-RG
<b>DESCRIPTION</b>	Core Polyurethane Chair – Medium Bench Height with Round Tube Base and Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21-28.5 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in



# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-PDHCH-MB-RG-T0-A0-RC
<b>DESCRIPTION</b>	Polyurethane Chair – Desk Height with Medium Back & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 16.25-21.5 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 15in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-PMBCH-MB-RG-T0-A0-CF-RC
<b>DESCRIPTION</b>	Polyurethane Chair – Medium Bench Height with Medium Back, Chrome Foot Ring, & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19.25-26.75 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 15in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-PHBCH-MB-RG-T0-A0-CF-RC
<b>DESCRIPTION</b>	Polyurethane Chair – High Bench Height with Medium Back & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 23.25-33.25 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 15in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS- PDHCH-MB-CR-T0-A0-CC
<b>DESCRIPTION</b>	Polyurethane Chair Chrome - Desk Height with Medium Back & Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17.5-22.75 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 15in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-PMBCH-MB-CR-T0-A0-CF-CC
<b>DESCRIPTION</b>	Polyurethane Chair Chrome - Medium Bench Height with Medium Back, Chrome Foot Ring, & Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.25-28 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 15in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-PHBCH-MB-CR-T0-A0-CF-CC
<b>DESCRIPTION</b>	Polyurethane Chair Chrome - High Bench Height with Medium Back, Chrome Foot Rin, & Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 24.5-34.75 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 15in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS- PDHCH-MB-RT-T0-A0-RG
<b>DESCRIPTION</b>	Polyurethane Chair with Round Tube Base - Desk Height with Medium Back & Stationary Glides in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19-24.25 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 15in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-BPMBCH-RT-RG
<b>DESCRIPTION</b>	Core Polyurethane Chair - Medium Bench Height with Round Tube Base and Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21-28.5 in
<b>SEAT</b>	Width: 16.5in, Depth: 15.5in, Thickness: 1.5in
<b>BACK</b>	Width: 14.5in, Height: 9in, Thickness: 1in



<b>CATALOG NO.</b>	MLS-HPDHCH-CR-T0-A0-CC
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair Chrome - Desk Height with Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17.75-23 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS- HPMBCH-CR-T0-A0-CF-CC
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair Chrome - Medium Bench Height with Chrome Foot Ring & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.5-28
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS- HPHBCH-CR-T0-A0-CF-CC
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair Chrome - High Bench Height with Chrome Foot Ring & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 24.75-35 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS- HPHBCH-CR-T0-A0-CF-CC
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair Chrome - High Bench Height with Chrome Foot Ring & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 24.75-35 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-HPDHCH-RT-T0-A0-RG
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair with Round Tube Base - Desk Height with Stationary Glides in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19.5-24.5 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-HPMBCH-RT-T0-A0-RG
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair with Round Tube Base - Medium Bench Height with Stationary Glides in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 22-29.5 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-HPDHCH-RG-T0-A0-RC
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair - Desk Height with Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17.25-22.5 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-HPMBCH-RG-T0-A0-CF-RC
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair - Medium Bench Height with Chrome Foot Ring & Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19.75-27.5 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-HPHBCH-RG-T0-A0-CF-RC
<b>DESCRIPTION</b>	Newport Industrial Polyurethane Chair - High Bench Height with Chrome Foot Ring & Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 24.25-34 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-TPDHCH-RG-A0-RC
<b>DESCRIPTION</b>	Erie Polyurethane Chair - Desk Height with Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 16.25-21 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-TPMBCH-RG-A0-CF-RC
<b>DESCRIPTION</b>	Erie Polyurethane Chair - Medium Bench Height with Chrome Foot Ring & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18.75-26.25 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-TPHBCH-RG-A0-CF-RC
<b>DESCRIPTION</b>	Erie Polyurethane Chair - High Bench Height with Chrome Foot Ring & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 23-33 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-TPDHCH-CR-A0-CC
<b>DESCRIPTION</b>	Erie Polyurethane Chair Chrome - Desk Height with Chrome Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 16.5 -21.75 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-TPMBCH-CR-A0-CF-CC
<b>DESCRIPTION</b>	Erie Polyurethane Chair Chrome – Medium Bench Height with Chrome Foot Ring & Chrome Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19.5–27 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-TPHBCH-CR-A0-CF-CC
<b>DESCRIPTION</b>	Erie Polyurethane Chair Chrome – High Bench Height with Chrome Foot Ring & Chrome Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 23.5–33.5 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-TPDHCH-RT-A0-RG
<b>DESCRIPTION</b>	Erie Polyurethane Chair – Desk Height with Round Tube Base & Stationary Glides in Gray Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18–23.25 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-TPMBCH-RT-A0-RG
<b>DESCRIPTION</b>	Erie Polyurethane Chair – Medium Bench Height with Round Tube Base & Stationary Glides in Gray Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21-28.5 in
<b>SEAT</b>	Width: 18in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 16.5in, Height: 12in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-PDHCH-RG-T0-A0-RC
<b>DESCRIPTION</b>	Polyurethane Chair – Desk Height with Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 16.25-21.5 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 10.25in, Thickness: 1 in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-PMBCH-RG-T0-A0-CF-RC
<b>DESCRIPTION</b>	Polyurethane Chair – Medium Bench Height with Chrome Foot Ring & Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19.25-26.75 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 10.25in, Thickness: 1 in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-PHBCH-RG-T0-A0-CF-RC
<b>DESCRIPTION</b>	Polyurethane Chair - High Bench Height with Chrome Foot Ring & Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 23.25-33.25 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 10.25in, Thickness: 1 in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-PDHCH-CR-T0-A0-CC
<b>DESCRIPTION</b>	Polyurethane Chair Chrome - Desk Height with Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17.5-22.75 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 10.25in, Thickness: 1 in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-PMBCH-CR-T0-A0-CF-CC
<b>DESCRIPTION</b>	Polyurethane Chair Chrome - Medium Bench Height with Chrome Foot Ring & Casters in Black Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.25-28 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 10.25in, Thickness: 1 in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY CHAIRS



<b>CATALOG NO.</b>	MLS-PHBCH-CR-T0-A0-CF-CC
<b>DESCRIPTION</b>	Polyurethane Chair Chrome - High Bench Height with Chrome Foot Ring & Casters in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 23.25-33.25 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 10.25in, Thickness: 1 in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-PDHCH-RT-T0-A0-RG
<b>DESCRIPTION</b>	Polyurethane Chair Round Tube Base - Desk Height with Stationary Glides in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19-24.25 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 10.25in, Thickness: 1 in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-PMBCH-RT-T0-A0-RG
<b>DESCRIPTION</b>	Polyurethane Chair Round Tube Base - Medium Bench Height with Stationary Glides in Blue Polyurethane
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21.5-29.25 in
<b>SEAT</b>	Width: 16in, Depth: 17in, Thickness: 1.5in
<b>BACK</b>	Width: 15.5in, Height: 10.25in, Thickness: 1 in
<b>COLOR</b>	Blue, Black, Grey

# LABORATORY STOOLS



<b>CATALOG NO.</b>	MLS- PTHBST-RG-CF-RG
<b>DESCRIPTION</b>	Polyurethane Tractor Sit-Stand with Chrome Foot Ring
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 22.5-32 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-PTDHSO-RG-RC
<b>DESCRIPTION</b>	Polyurethane Tractor Sit-Stand Stool - Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 16.5-21.5 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-PTHBST-CR-CF-RG
<b>DESCRIPTION</b>	Polyurethane Tractor Sit-Stand Chrome with Chrome Foot Ring
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 23-33.5
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS- PTMBSO-RG-CF-RC
<b>DESCRIPTION</b>	Polyurethane Tractor Sit-Stand with Chrome Foot Ring
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19-26.5 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey



# LABORATORY STOOLS



<b>CATALOG NO.</b>	MLS-PTDHSO-CR-CC
<b>DESCRIPTION</b>	Polyurethane Tractor Sit-Stand Stool Chrome - Desk Height with Chrome Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17.5-22.5 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-PTMBSO-CR-CF-CC
<b>DESCRIPTION</b>	Polyurethane Tractor Sit-Stand Stool Chrome - Medium Bench Height with Chrome Foot Ring & Chrome Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20-28 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-PTDHSL-RG-RC
<b>DESCRIPTION</b>	Polyurethane Tractor Stool - Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 16-21 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS- PTMBSL-RG-CF-RC
<b>DESCRIPTION</b>	Polyurethane Tractor Stool - Medium Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19-26.5 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey





# LABORATORY STOOLS



<b>CATALOG NO.</b>	MLS-PTDHSL-CR-CC
<b>DESCRIPTION</b>	Polyurethane Tractor Stool Chrome – Desk Height with Chrome Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17-22 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-PTMBSL-CR-CF-CC
<b>DESCRIPTION</b>	Polyurethane Tractor Stool Chrome – Medium Bench Height with Chrome Foot Ring & Chrome Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20-27.5 in
<b>SEAT</b>	Width: 17in, Depth: 14in, Thickness: 2.5in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-IPDHSO-RT-RG
<b>DESCRIPTION</b>	Huron Polyurethane Stool with Round Tube Base – Desk Height with Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18-23 in
<b>SEAT</b>	Width: 13in, Depth: 13in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-IPMBSO-RT-RG
<b>DESCRIPTION</b>	Huron Polyurethane Stool with Round Tube Base – Medium Bench Height with Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.5-28.5
<b>SEAT</b>	Width: 13in, Depth: 13in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



# LABORATORY STOOLS



<b>CATALOG NO.</b>	MLS-IPDHSO-RG-RC
<b>DESCRIPTION</b>	Huron Polyurethane Stool - Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 15.5-20.5 in
<b>SEAT</b>	Width: 13in, Depth: 13in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-IPMBSO-RG-CF-RC
<b>DESCRIPTION</b>	Huron Polyurethane Stool - Medium Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17.75-25.5 in
<b>SEAT</b>	Width: 13in, Depth: 13in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey



<b>CATALOG NO.</b>	MLS-IPHBSO-RG-CF-RC
<b>DESCRIPTION</b>	Huron Polyurethane Stool - High Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 22-32 in
<b>SEAT</b>	Width: 13in, Depth: 13in, Thickness: 1in
<b>COLOR</b>	Blue, Black, Grey

<b>CATALOG NO.</b>	MLS-VDHSO-RG-RC-8613
<b>DESCRIPTION</b>	Vinyl Stool without Back - Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 17.25-22.5 in
<b>SEAT</b>	Diameter 17in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



# LABORATORY STOOLS



<b>CATALOG NO.</b>	MLS-VMBSO-RG-CF-RC-8613
<b>DESCRIPTION</b>	Vinyl Stool without Back - Medium Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19.75-27.75 in
<b>SEAT</b>	Diameter 17in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

<b>CATALOG NO.</b>	MLS-VHBSO-RG-CF-RC-8613
<b>DESCRIPTION</b>	Vinyl Stool without Back - High Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 24.5-34.5 in
<b>SEAT</b>	Diameter 17in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



<b>CATALOG NO.</b>	MLS-VDHSO-CR-CC-8801
<b>DESCRIPTION</b>	Vinyl Stool without Back Chrome - Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 18.25-23.25 in
<b>SEAT</b>	Diameter 17in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

# LABORATORY STOOLS



<b>CATALOG NO.</b>	MLS-VMBSO-CR-CF-CC-8801
<b>DESCRIPTION</b>	Vinyl Stool without Back Chrome - Medium Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21.25-28.75 in
<b>SEAT</b>	Diameter 17in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

<b>CATALOG NO.</b>	MLS-VHBSO-CR-CF-CC-8801
<b>DESCRIPTION</b>	Vinyl Stool without Back Chrome - High Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 25.5-35.25 in
<b>SEAT</b>	Diameter 17in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.



<b>CATALOG NO.</b>	MLS-VDHSO-RT-RG-8613
<b>DESCRIPTION</b>	Vinyl Stool without Back - Desk Height with Round Tube Base & Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 19.75-25.25 in
<b>SEAT</b>	Diameter 17in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

# LABORATORY STOOLS



<b>CATALOG NO.</b>	MLS-VMBSO-RT-RG-8613
<b>DESCRIPTION</b>	Vinyl Stool without Back - Medium Bench Height with Round Tube Base & Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 22.5-30.25
<b>SEAT</b>	Diameter 17in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

<b>CATALOG NO.</b>	MLS-WDHSO-RG-RC
<b>DESCRIPTION</b>	Wood Stool - Desk Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 16-21 in
<b>SEAT</b>	Diameter 13.75in, Thickness: 0.75in



<b>CATALOG NO.</b>	MLS-WHBSO-RG-CF-RC
<b>DESCRIPTION</b>	Wood Stool - High Bench Height with Chrome Foot Ring & Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 22.25-32 in
<b>SEAT</b>	Diameter 13.75in, Thickness: 0.75in

<b>CATALOG NO.</b>	MLS-WMBSO-RT-RG
<b>DESCRIPTION</b>	Wood Stool - Medium Bench Height with Round Tube Base & Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.5-28 in
<b>SEAT</b>	Diameter 13.75in, Thickness: 0.75in





# LABORATORY STOOLS



<b>CATALOG NO.</b>	MLS-VSLMD-CR-RC-8801
<b>DESCRIPTION</b>	Vinyl Screw-Lift Stool – Medium Bench Height with Casters
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.5-28 in
<b>SEAT</b>	Diameter 14.5in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

<b>CATALOG NO.</b>	MLS-PSTSS-RG-CF-RG
<b>DESCRIPTION</b>	Polyurethane Sit-Stand – Standing Height with Chrome Foot Ring
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 21.5-31.5 in
<b>SEAT</b>	Width: 13in, Depth: 10.5in, Thickness: 2.5in



<b>CATALOG NO.</b>	MLS-PSTSS-CR-CF-RG
<b>DESCRIPTION</b>	Polyurethane Sit-Stand Chrome – Standing Height with Chrome Foot Ring
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 22.5-32.5 in
<b>SEAT</b>	Width: 13in, Depth: 10.5in, Thickness: 2.5in

<b>CATALOG NO.</b>	MLS-VMBMS-RG-NF-RG-8613
<b>DESCRIPTION</b>	Vinyl Mini-Stool – Medium Bench Height with Stationary Glides
<b>ADJUSTMENT</b>	Pneumatic Adjustment: 20.75-27.75 in
<b>SEAT</b>	Diameter 14.5in, Thickness: 3.5in
<b>COLOR</b>	Orange Kist, Black, Blue, Burgundy, Charcoal, Taupe, Imperial, Dove, Adobe, Storm, Taupe, Borscht, Carbon, Sterling, Marine.

