

# Laboratory Plasticware



# CENTRIFUGE TUBES



CATALOG NO.	DESCRIPTION
MCT-00015	15 ml Conical Bottom, Grad, Flat Cap, EO
MCT-00050	50 ml Conical Bottom, Grad, Flat Cap, EO
<b>Polypropylene tubes include blue HDPE leak-proof screw caps. Highly smooth hydrophobic internal surfaces minimize wall effect. Silk-screen graduations and a frosted surface make writing easier.</b>	

## Features:

- High strength non-toxic materials
- Raw material comply with USP VI requirements
- Printed clear graduations with marking area
- One-hand operation design
- DNase, RNase and pyrogen free

CATALOG NO.	DESCRIPTION
MLS-3131-345-008	Labcon 15mL SuperClear Centrifuge Tubes
MLS-3181-345-008	Labcon 50 mL Super Clear Centrifuge Tubes
<b>SuperClear® Screw Cap Micro Tubes with molded-in elastomeric seal are ideal for any sample processing application. Excellent for frozen storage or when working with hazardous or valuable materials.</b>	

## Features:

- Elastomeric polyethylene caps that prevent leaking - Freeze to -80 °C
- Clear graduations and resin allow for accurate sample observation
- Non-Sterile or Sterile; Bulk packaging (50/bag, 10 bags/case, 500 tubes/case)
- Biologically Tested to be certified free of RNase/DNase, DNA, ATP, Endotoxins, and Pyrogenic Materials



CATALOG NO.	DESCRIPTION
MLS-D1002	Self Standing centrifuge tubes, PP/HDPE, 50 ml, Starline
<b>Polypropylene tubes feature a convenient flat bottom so they stand on their own. Include blue HDPE leak-proof screw caps, silk-screen black graduations and a large white marking spot.</b>	

# CENTRIFUGE TUBES



CATALOG NO.	DESCRIPTION
MLS-3095-340-008	50mL PerformR® Freestanding Centrifuge Tubes
<b>PerformR® freestanding centrifuge tubes are designed for your everyday centrifugation and storage needs. These tubes include large white solvent resistant labeling areas and black easy to read graduations.</b>	

CATALOG NO.	DESCRIPTION
MLS-3144-335-001	Labcon 15 mL Centrifuge Tube
MLS-3184-345-001	Labcon 50 mL Centrifuge Tube
<b>The MetalFree® tubes are made with ultra clean resins for the tube and cap and processed to ensure that common trace metals are kept below an ICP-MS detection level of one part per billion (1ppb). These tubes are perfect for storage of samples and for applications requiring the assurance of low background metals.</b>	



CATALOG NO.	DESCRIPTION
MLS-10482	15 ml Centrifuge tubes, Sterile
MLS-10485	50 ml Centrifuge tubes, Sterile
MLS-10480	15 ml Centrifuge tubes, Non-Sterile
MLS-10483	50 ml Centrifuge tubes, Non-Sterile
<b>Manufactured using high-purity, virgin USP Class VI Medical grade polypropylene, these centrifuge tubes are free from natural rubber and heavy metal and conform to US FDA 21 CFR. Include leak-proof threaded HDPE caps that are free of added colorings and trace metals.</b>	

# MICRO CENTRIFUGE TUBES



CATALOG NO.	DESCRIPTION
MLS-MCT1.5	1.5 ml Micro centrifuge tubes.
<b>These Micro Centrifuge Tubes are Conical, Graduated and Non-sterilized. These come with attached caps.</b>	

CATALOG NO.	DESCRIPTION
MLS-111642	Microcentrifuge Tube, 1.5mL, PP, No Cap, Graduated, Natural
MLS-111620	Microcentrifuge Tube, 2mL, PP, No Cap, Graduated, Natural
<b>These polypropylene (PP) microcentrifuge tubes without cap are resistant to most chemical reagents and can be autoclaved at 130°C.</b>	



CATALOG NO.	DESCRIPTION
MLS-111550	600uL, Microcentrifuge Tube, PP, with Attached Click-Seal Snap Cap, Natural
MLS-110550	500uL, Microcentrifuge Tube, PP, with Attached Click-Seal Snap Cap, Natural
<b>These polypropylene microcentrifuge tubes feature attached leak-resistant flat snap caps. Each tube is molded with precise graduation marks for accurate measurements.</b>	

CATALOG NO.	DESCRIPTION
MLS-111568	2mL Microcentrifuge Tube, PP, Attached Snap Cap, Graduated, Natural
<b>These polypropylene microcentrifuge tubes feature attached leak-resistant flat snap caps. Each tube is molded with precise graduation marks for accurate measurements.</b>	





# MICRO CENTRIFUGE TUBES



CATALOG NO.	DESCRIPTION
MLS-111552LK	Microcentrifuge Tube, 0.5mL, PP, Attached Locking Snap Cap.
MLS-111562LK	Microcentrifuge Tube, 1.5mL, PP, Attached Locking Snap Cap, Graduated, Natural
MLS-111572LK	Microcentrifuge Tube, 2.0mL, PP, Attached Locking Snap Cap, Graduated, Natural
MLS-111582LK	Microcentrifuge Tube, 5.0mL, PP, Attached Locking Snap Cap, Graduated, Natural
<b>These tubes are engineered with a secure locking system that provides the security of a screw cap microtube with the simplicity of a snap cap tube.</b>	

CATALOG NO.	DESCRIPTION
MLS-110540	500uL (0.5mL) Microcentrifuge Tube, 6 x 44mm, PE, Attached Plug Cap
<b>These 500uL (0.5mL) Beckman® Type microcentrifuge tubes with attached cap are non-graduated. The tubes are made from polyethylene and are not autoclavable.</b>	



CATALOG NO.	DESCRIPTION
MLS-111554A	Microcentrifuge Tube, 0.5mL, PP, Attached Snap Cap, Graduated, AMBER
MLS-111564A	1.5mL Microcentrifuge Tube, PP, Attached Snap Cap, AMBER
MLS-111574A	Microcentrifuge Tube, 2.0mL, PP, Attached Snap Cap, Graduated, AMBER
<b>Manufactured in a class 100,000 cleanroom and are certified to be free from RNase, DNase, DNA, and Pyrogens</b>	



# MICRO CENTRIFUGE TUBES



CATALOG NO.	DESCRIPTION
MLS-111554B	Microcentrifuge Tube, 0.5mL, PP, Attached Snap Cap, Graduated, BLUE
MLS-111554G	Microcentrifuge Tube, 0.5mL, PP, Attached Snap Cap, Graduated, GREEN
MLS-111554R	Microcentrifuge Tube, 0.5mL, PP, Attached Snap Cap, Graduated, RED
MLS-111554Y	Microcentrifuge Tube, 0.5mL, PP, Attached Snap Cap, Graduated, YELLOW

CATALOG NO.	DESCRIPTION
MLS-111564B	1.5mL Microcentrifuge Tube, PP, Attached Snap Cap, BLUE
MLS-111564G	1.5mL Microcentrifuge Tube, PP, Attached Snap Cap, GREEN
MLS-111564R	1.5mL Microcentrifuge Tube, PP, Attached Snap Cap, RED
MLS-111564Y	1.5mL Microcentrifuge Tube, PP, Attached Snap Cap, YELLOW

CATALOG NO.	DESCRIPTION
MLS-111574B	Microcentrifuge Tube, 2.0mL, PP, Attached Snap Cap, Graduated, BLUE
MLS-111574G	Microcentrifuge Tube, 2.0mL, PP, Attached Snap Cap, Graduated, GREEN
MLS-111574R	Microcentrifuge Tube, 2.0mL, PP, Attached Snap Cap, Graduated, RED
MLS-111574Y	Microcentrifuge Tube, 2.0mL, PP, Attached Snap Cap, Graduated, YELLOW

# PETRI DISHES



CATALOG NO.	DESCRIPTION
MPD-00060-500	Clear PS, 60x15mm, Sterile
MPD-00100-500	Clear PS, 100x15mm, Sterile
MPD-00150-500	Clear PS, 150x15mm, Sterile
<b>These petri dishes are clear made up of polystyrene, 3 vents.</b>	

## Features:

- Flat transparent surface
- Numeric mark on the bottom for easy identification
- PA/PE packaging
- Vent for gas exchange
- Stackable for easy storage and handling
- Sterilized by E-Beam, SAL=10<sup>-6</sup>
- Non-pyrogenic

# SAMPLE CONTAINER

An IVD CERTIFIED PRODUCTS best recommended for diagnostics with a leakproof screw cap, principally important when transporting hazardous material. Samples containers are shatter resistant, Manufactured from ultra-pure USP-VI Virgin Polypropylene (PP) and Comply with the requirement of 21 CFR 177.1520

CATALOG NO.	DESCRIPTION
MLS-SC-120	120 ml, Individually wrapped, Sterile, Frosted portion for easy writing, Leak proof screw cap
MLS-C-11008-N	8 oz, polypropylene
<b>Sterile (Gamma Irradiation), DNase RNase Free, Human DNA Free, Free of heavy metals, Leak-proof Tested 635 mm/Hg, ATP &amp; Pyrogen free, Autoclavable</b>	



## Features:

- Ultra clear polypropylene container.
- Frosted portion for easy writing.
- Clear and easy to read graduations.
- Sterile container is individually packed in polyethylene pouch.

# SEROLOGICAL PIPETS



CATALOG NO.	DESCRIPTION
MSP-00001	1ml Serological Pipets
MSP-00002	2ml Serological Pipets
MSP-00005	5ml Serological Pipets
MSP-00010	10ml Serological Pipets
MSP-00025	25ml Serological Pipets
MSP-00050	50ml Serological Pipets with black stripe
<b>These Serological Pipets are Individually Wrapped and Sterile</b>	

## Features:

- *Edged color-coded band for identification*
- *Negative graduations for extra volume*
- *Sterilized by E-beam*
- *Non-pyrogenic, non-cytotoxic and non-haemolytic*
- *DNase/RNase free*
- *Polyolefin fiber filter to reduce contamination*

# PIPETTE FILTER TIPS

CATALOG NO.	DESCRIPTION
MLS-FPT-0010	10ul Filter Tips, graduated, Natural color, Sterile, Racked
MLS-FPT-0020	20ul Filter Tips, graduated, Natural color, Sterile, Racked
MLS-FPT-0200	200ul Filter Tips, graduated, Natural color, Sterile, Racked
MLS-FPT-1000	1000ul Filter Tips, graduated, Natural color, Sterile, Racked
MLS-FPT-1250	1250ul Filter Tips, graduated, Natural color, Sterile, Racked
<b>These low Retention, Racked are disposable tips intended for use with most brands of pipettors. Manufactured from 100% pure porous polyethylene, filters are designed to prevent liquids and aerosols from entering the body of the micropipette to protect against cross-contamination.</b>	





# PIPETTE TIPS

Made from medical grade PP and manufactured in ISO 8 grade cleaning room. With precision mould and professional injection skill, these tips has smooth inner surface which is very important to ensure accurate results during liquid handling.

CATALOG NO.	DESCRIPTION
MPT-00010N	10uL(0.1uL-10uL), Pipette Tips, Bulk Non Sterile
MPT-00200N	200uL(10uL-200uL), Pipette Tips, Bulk Non Sterile
MPT-01000N	1000uL(100uL-1000uL), Pipette Tips, Bulk Non Sterile
MLS-PT1250	1250 ul Pipette Tip, Natural, Sterile.



## Racked Pipette Tips

Best suited for biological applications, these tips are non-cytotoxic, non-hemolytic, and have 99.9% Bacterial Aerosol Filtration Efficiency (BFE). Sterile tips are manufactured from ultrapure USP Class VI certified, virgin polypropylene. DNase, RNase and Pyrogen free.



CATALOG NO.	DESCRIPTION
MLS-RPT-0010	10ul Pipette Tips, graduated, Natural color, Sterile, Racked
MLS-RPT-0200	200ul Pipette Tips, graduated, Natural color, Sterile, Racked
MLS-RPT-1000	1000ul Pipette Tips, graduated, Natural color, Sterile, Racked
MLS-RPT-1250	1250ul Pipette Tips, graduated, Natural color, Sterile, Racked

# STORAGE BOXES

CATALOG NO.	OUTSIDE DIMENSIONS	INSIDE DIMENSIONS	CAPACITY
MLS-S-14599	18 x 12 x 12"	14 x 9 x 10"	30
MLS-S-14600	26 x 16 x 14"	20 x 12 x 11"	70
MLS-S-14601	33 x 20 x 14"	28 x 17 x 11"	116
Comfortable, innovative recessed latch for easy handling. Stackable when storing product. Durable polypropylene construction.			



# TRANSFER PIPETS



CATALOG NO.	MLS-134010	MLS-134020	MLS-134070	MLS-134050	MLS-134060	MLS-134090
CAPACITY	1.3mL	1.5mL	5mL	5.8mL	5.8mL	8.7mL
DROPS	45mL	50mL	50mL	60mL	50mL	45mL
BULB DRAW	1mL	1mL	3.2mL	3.4mL	3.4mL	4.8mL
TOTAL LENGTH	51mm	104mm	153mm	147mm	157mm	147mm
STEM DIAMETER	1.15mm	4.5mm	4.6mm	4.6mm	4.6mm	4.6mm
TIP DIAMETER	1.15mm	0.9mm	0.85mm	1.00mm	0.85mm	1.2mm
FINE TIP TRANSFER PIPETS ARE IDEAL FOR DISPENSING SMALL VOLUMES, LOADING HEMOCYTOMETERS AND TRANSFERRING LIQUID TO AND FROM MULTI-WELL PLATES.						



CATALOG NO.	DESCRIPTION
MLS-P31205	Plastic transfer pipettes, LDPE 3mL, Graduated, 154mm height
DNase RNase Free, Free of endotoxins, Free of heavy metals, ATP & Pyrogen free, Human DNA free, Sterile (Gamma Irridiation).	

# TRANSFER PIPETS



CATALOG NO.	MLS-138020	MLS-138060	MLS-138050	MLS-138080	MLS-138030	MLS-138090
CAPACITY	1.7mL	4mL	5mL	7mL	8mL	8.5mL
DROPS	25 ml	25 ml	23 ml	20 ml	20 ml	18 ml
BULB DRAW	0.9mL	2mL	2mL	3.2mL	4.6mL	4.9mL
TOTAL LENGTH	87mm	130mm	155mm	155mm	152mm	137mm
STEM DIAMETER	4.5mm	6.6mm	6.8mm	7.8mm	6.2mm	6.1mm
TIP DIAMETER	2.4mm	2.5mm	2.4mm	3.1mm	3mm	4mm
OUR NON-GRADUATED TRANSFER PIPETS ARE OFFERED IN A VARIETY OF SIZES AND CAPACITIES						



CATALOG NO.	MLS-138040	MLS-136030	MLS-139020	MLS-139040
CAPACITY	1.2mL	4mL	6.5mL	15mL
DROPS	25mL	24mL	28mL	19mL
BULB DRAW	0.9mL	3mL	3.2mL	8mL
TOTAL LENGTH	65mm	85mm	155mm	155mm
STEM DIAMETER	2.5mm	2.5mm	7.8mm	3mm
TIP DIAMETER	2.5mm	2.5mm	2.1mm	3mm
THESE TRANSFER PIPETS FEATURE LONG, THIN STEMS THAT ARE FLEXIBLE AND CAN BE BENT TO ACCESS NARROW OPENINGS AS WELL AS SMALL DIAMETER TUBES.				

# TRANSFER PIPETS



CATALOG NO.	MLS-136036	MLS-137035	MLS-137010	MLS-138070	MLS-137050	MLS-135030
CAPACITY	1.5mL	3mL	5mL	5mL	5mL	7mL
DROPS	21 ml	25 ml	25 ml	21 ml	23 ml	20 ml
BULB DRAW	0.9mL	2.2mL	3.4mL	1.8mL	2mL	3.2mL
TOTAL LENGTH	116mm	140mm	148mm	150mm	155mm	160mm
STEM DIAMETER	3.2mm	4.4mm	4.8mm	7.2mm	6.8mm	7.9mm
TIP DIAMETER	3.2mm	2.4mm	2.7mm	3mm	2.6mm	3mm
GRADUATED TO	0.3mL	1mL	1mL	2.25mL	2mL	3mL
GRADUATED TRANSFER PIPETS ARE DESIGNED FOR MEASUREMENTS, THESE PIPETS FEATURE MOLDED-IN GRADUATION MARKS THAT WILL NEVER FADE OR WASH OFF.						

## INOCULATING LOOPS



CATALOG NO.	DESCRIPTION
MIL-00010	10 uL Blue, Sterile, Soft
MIL-00001	1 uL White, Sterile, Soft

### Features:

- Polypropylene material
- Hydrophilic treated
- Flexible handles for easier collection and inoculation
- Sterilized by Ethylene Oxide

# SYRINGE FILTERS



## Nylon Syringe Filters

CATALOG NO.	DESCRIPTION	PORE SIZE	HOUSING DIAMETER
MSF-NY2213	Nonsterile Nylon Syringe Filters	0.22 µm	13 mm
MSF-NY4513	Nonsterile Nylon Syringe Filters	0.45 µm	13 mm
MSF-NY2225	Nonsterile Nylon Syringe Filters	0.22 µm	25 mm
MSF-NY4525	Nonsterile Nylon Syringe Filters	0.45 µm	25 mm

## Hydrophobic PTFE Syringe Filters

CATALOG NO.	DESCRIPTION	PORE SIZE	HOUSING DIAMETER
MSF-PT2213	Nonsterile Hydrophobic PTFE Syringe Filters	0.22 µm	13 mm
MSF-PT4513	Nonsterile Hydrophobic PTFE Syringe Filters	0.45 µm	13 mm
MSF-PT2225	Nonsterile Hydrophobic PTFE Syringe Filters	0.22 µm	25 mm
MSF-PT4525	Nonsterile Hydrophobic PTFE Syringe Filters	0.45 µm	25 mm

## Hydrophobic PVDF Syringe Filters

CATALOG NO.	DESCRIPTION	PORE SIZE	HOUSING DIAMETER
MSF-PV2213	Nonsterile Hydrophobic PVDF Syringe Filters	0.22 µm	13 mm
MSF-PV4517	Nonsterile Hydrophobic PVDF Syringe Filters	0.45 µm	13 mm
MSF-PV2225	Nonsterile Hydrophobic PVDF Syringe Filters	0.22 µm	25 mm
MSF-PV4525	Nonsterile Hydrophobic PVDF Syringe Filters	0.45 µm	25 mm



# SYRINGE FILTERS

## PES Syringe Filters



CATALOG NO.	DESCRIPTION
MSF-PES-0825	PES Syringe Filter, Sterile, 8.0 Micron, 25mm
MSF-PES-4528	PES Syringe Filters, Sterile, 0.45 Micron, 28mm
MSF-PES-2028	PES Syringe Filters, Sterile, 0.20 Micron, 28mm
MSF-PES-2213	PES Syringe Filters, 0.22 Micron, 13mm
MSF-PES-4513	PES Syringe Filters, 0.45 Micron, 13mm
MSF-PES-2217	PES Syringe Filters, 0.22 Micron, 17mm,

## Vent Filters With PTFE Membranes



CATALOG NO.	DESCRIPTION
MSF-VF-02713	PTFE Disposable Vent Filter, 0.2 Micron, 50mm, 7-13mm Hose Barb Connection
MSF-VF-02612	PTFE Vent Filter, Sterile, 0.2 Micron, 50mm, 6 To 12mm Hose BARB
MSF-VF-04550	PTFE Vent Filter, Sterile, 0.45 Micron, 50mm, 1/4 In To 3/8 In Hose BARB

# BEAKERS

## Plastic Beakers

Polymethylpentene beakers with pouring spouts are ideal for general laboratory use. Tapered design is convenient for stacking and holding. All sizes have easy to read silk-screened graduations.



CATALOG NO.	DESCRIPTION
MLS-PMP-50	50 ml, Beakers, Printed graduations, PMP
MLS-PMP-100	100 ml, Beakers, Printed graduation, PMP
MLS-PMP-250	250 ml, Beakers, Printed graduation, PMP
MLS-PMP-500	500 ml, Beakers, Printed graduation, PMP
MLS-PMP-1000	1000 ml, Beakers, Printed graduation, PMP
MLS-PMP-2000	2000 ml, Beakers, Printed graduation, PMP
MLS-PMP-5000	5000 ml, Beakers, Printed graduation, PMP
MLS-PMP-10000	10000 ml, Beakers, Printed graduation, PMP

## Beakers with Handle

Polymethylpentene beakers with pouring spouts are ideal for general laboratory use. Includes rectangular handle for easy carrying and pouring. All sizes have easy to read silk-screened graduations.

CATALOG NO.	DESCRIPTION
MLS-PBH-100	100 ml, Beakers with Handle, Printed graduation, PMP
MLS-PBH-250	250 ml, Beakers with Handle, Printed graduation, PMP
MLS-PBH-500	500 ml, Beakers with Handle, Printed graduation, PMP
MLS-PBH-1000	1000 ml, Beakers with Handle, Printed graduation, PMP
MLS-PBH-2000	2000 ml, Beakers with Handle, Printed graduation, PMP
MLS-PBH-5000	5000 ml, Beakers with Handle, Printed graduation, PMP
MLS-PBH-10000	10000 ml, Beakers with Handle, Printed graduation, PMP



# BEAKERS

## Wide Spout Beakers

Beakers, molded in polypropylene, have excellent clarity and very good chemical resistance. They are autoclavable and feature raised graduations on the outer surface. Beakers feature oversized, tapered spouts for easier pouring of solutions.



CATALOG NO.	DESCRIPTION
MLS-PP-0025	Beakers, Wide spout, PP, 25 ml
MLS-PP-0050	Beakers, Wide spout, PP, 50 ml
MLS-PP-0100	Beakers, Wide spout, PP, 100 ml
MLS-PP-0250	Beakers, Wide spout, PP, 250 ml
MLS-PP-0500	Beakers, Wide spout, PP, 500 ml
MLS-PP-1000	Beakers, Wide spout, PP, 1000 ml
MLS-PP-2000	Beakers, Wide spout, PP, 2000 ml
MLS-PP-5000	Beakers, Wide spout, PP, 5000 ml

## Griffin Style Beakers

Beakers molded in polypropylene have excellent clarity and chemical resistance. They are autoclavable and feature raised graduations on the outer surface. Griffin Style beakers nest without sticking together.

CATALOG NO.	DESCRIPTION
MLS-PPG-0050	Beakers, Griffin Style, Polypropylene, 50 ml
MLS-PPG-0100	Beakers, Griffin Style, Polypropylene, 100 ml
MLS-PPG-0250	Beakers, Griffin Style, Polypropylene, 250 ml
MLS-PPG-0500	Beakers, Griffin Style, Polypropylene, 500 ml
MLS-PPG-0600	Beakers, Griffin Style, Polypropylene, 600 ml
MLS-PPG-1000	Beakers, Griffin Style, Polypropylene, 1000 ml
MLS-PPG-2000	Beakers, Griffin Style, Polypropylene, 2000 ml



# BEAKERS

## Tri-Corner Beakers

These tri-pour plastic beakers feature dripless pouring spouts and range in size from 50mL to 1000mL. They are unbreakable and can be used with commonly used acids, solvents and alkalis. The laboratory beakers can be autoclaved and are reusable, yet are inexpensive enough to be disposable. All six sizes are graduated in mL.



CATALOG NO.	DESCRIPTION
MLS-PPT-0050	Beaker, Three Corner, PP, Graduated, 50mL
MLS-PPT-0100	Beaker, Three Corner, PP, Graduated, 100mL
MLS-PPT-0250	Beaker, Three Corner, PP, Graduated, 250mL
MLS-PPT-0400	Beaker, Three Corner, PP, Graduated, 400mL
MLS-PPT-0800	Beaker, Three Corner, PP, Graduated, 800mL
MLS-PPT-1000	Beaker, Three Corner, PP, Graduated, 1000mL


## Tri-Corner Beakers

These tri-pour plastic beakers range from 50mL to 1000mL. They are unbreakable and can be used with commonly used acids, solvents and alkalis. The laboratory beakers can be autoclaved and are reusable, yet are inexpensive enough to be disposable. All six sizes are graduated in mL.

CATALOG NO.	DESCRIPTION
MLS-CPP-0050	Beaker, Three Corner, Clarified PP, Graduated, 50mL
MLS-CPP-0100	Beaker, Three Corner, Clarified PP, Graduated, 100mL
MLS-CPP-0250	Beaker, Three Corner, Clarified PP, Graduated, 250mL
MLS-CPP-0400	Beaker, Three Corner, Clarified PP, Graduated, 400mL
MLS-CPP-0800	Beaker, Three Corner, Clarified PP, Graduated, 800mL
MLS-CPP-1000	Beaker, Three Corner, Clarified PP, Graduated, 1000mL



# MEASURING CYLINDERS



CATALOG NO.	DESCRIPTION
MLS-CPSET7-H	Graduated Cylinder Set of 7, PP, Hex-Base
MLS-CPSET5	Graduated Cylinder Set of 5, PP, Hex-Base
<b>A Convenient set of 7 graduated polypropylene cylinders in sizes 10ml, 25ml, 50ml, 100ml, 250ml, 500ml, and 1000ml. Cylinders feature hexagonal base that provides excellent stability. Also available as a set of five in sizes 10ml, 50ml, 100ml, 250ml and 500ml.</b>	



CATALOG NO.	DESCRIPTION
MLS-81150	Conical Measuring Cylinders, PP, 12 ML
MLS-81151	Conical Measuring Cylinders, PP, 25 ML
MLS-81152	Conical Measuring Cylinders, PP, 50mL
MLS-81153	Conical Measuring Cylinders, PP, 125mL
MLS-81154	Conical Measuring Cylinders, PP, 200mL
<b>Our polypropylene conical measuring cylinders are ideal for measuring small quantities of reagents. They are autoclavable and have excellent contact clarity. The hexagonal base provides excellent stability to the cylinder, and the tapered spout enables easy pouring of reagents. Features easy-to-read raised graduations.</b>	



# MEASURING CYLINDERS



CATALOG NO.	DESCRIPTION
MLS-P50501	Measuring Cylinder, PMP, Class A,, 10 ml
MLS-P50502	Measuring Cylinder, PMP, Class A, 25 ml
MLS-P50503	Measuring Cylinder, PMP, Class A, 50 ml
MLS-P50504	Measuring Cylinder, PMP, Class A, 100 ml
MLS-P50505	Measuring Cylinder, PMP, Class A, 250 ml
MLS-P50506	Measuring Cylinder, PMP, Class A, 500 ml
MLS-P50507	Measuring Cylinder, PMP, Class A, 1000 ml
MLS-P50508	Measuring Cylinder, PMP, Class A, 2000 ml
<b>Meets accuracy requirements of Class A under DIN 12681. Supplied with conformity test certificate that indicates batch number and year of production.</b>	

CATALOG NO.	DESCRIPTION
MLS-P50401	Measuring Cylinder, PMP, Class B, 10 ml
MLS-P50402	Measuring Cylinder, PMP, Class B, 25 ml
MLS-P50403	Measuring Cylinder, PMP, Class B, 50 ml
MLS-P50404	Measuring Cylinder, PMP, Class B, 100 ml
MLS-P50405	Measuring Cylinder, PMP, Class B, 250 ml
MLS-P50406	Measuring Cylinder, PMP, Class B, 500 ml
MLS-P50407	Measuring Cylinder, PMP, Class B, 1000 ml
MLS-P50408	Measuring Cylinder, PMP, Class B, 2000 ml
<b>Meets accuracy requirements of Class A under DIN 12681. Supplied with conformity test certificate that indicates batch number and year of production.</b>	



# MEASURING CYLINDERS



CATALOG NO.	DESCRIPTION
MLS-80031	Measuring Cylinder, Hex. Base, PP, 10 ml
MLS-80032	Measuring Cylinder, Hex. Base, PP, 25 ml
MLS-80033	Measuring Cylinder, Hex. Base, PP, 50 ml
MLS-80034	Measuring Cylinder, Hex. Base, PP, 100 ml
MLS-80035	Measuring Cylinder, Hex. Base, PP, 250 ml
MLS-80036	Measuring Cylinder, Hex. Base, PP, 500 ml
MLS-80037	Measuring Cylinder, Hex. Base, PP, 1000 ml
MLS-80038	Measuring Cylinder, Hex. Base, PP, 2000 ml
<b>Meets accuracy requirements of Class A under DIN 12681. Supplied with conformity test certificate that indicates batch number and year of production.</b>	

## ACRYLIC BURET



CATALOG NO.	DESCRIPTION
MLS-11501	Acrylic Buret, 25 ml
MLS-11502	Acrylic Buret, 50 ml
MLS-11503	Acrylic Buret, 100 ml
<b>This buret has a clear acrylic body and leak-proof stopcock with a selflubricating PTFE plug. Acrylic body is resistant to dilute mineral acids and bases except Hydrofluoric Acid and Ammonium Hydroxide. Graduations are in 0.1ml increments.</b>	

# BOTTLES

## Dispenser Dropping Bottles



CATALOG NO.	DESCRIPTION
MLS-DDB015	Dispenser dropping bottle 15mL
MLS-DDB030	Dispenser dropping bottle 30mL
MLS-DDB060	Dispenser dropping bottle 60mL
MLS-DDB125	Dispenser dropping bottle 125mL
MLS-DDB250	Dispenser dropping bottle 250mL
<b>This LDPE dispensing bottle can be easily operated with one hand. The bottle features a polypropylene screw cap with a tethered spout lid that snaps securely in place.</b>	

## Aspirator Bottles

CATALOG NO.	DESCRIPTION
MLS-34101	Aspirator Bottles, PP 5-Liter
MLS-34102	Aspirator Bottles, PP 10-Liter
MLS-34103	Aspirator Bottles, PP 20-Liter
<b>The leak-proof spigot screwed to the bottle requires only a half turn to open or close. Polyethylene cap and stopcock should be removed before autoclaving.</b>	



## Drosophila Bottles



CATALOG NO.	DESCRIPTION
MLS-DB-6OZ	Drosophila Bottles
<b>Made of high quality polypropylene with good visibility. The bottle capacity is 6 oz (177ml) with dimensions of 57L x 57W x 103mmH. Square Bottom, Disposable, and Autoclavable. Graduations are molded.</b>	

# BOTTLES

## Wash Bottles



These wash bottles are made from industrial grade low-density polyethylene with polypropylene stem closure and venting assembly. They are pre-printed with the Chemical name and formula, the National Fire Protection Association (NFPA) Section 704 Four-color Diamond, health hazards, Chemical Abstract System (CAS) number, and suggested protective clothing and equipment. Each bottle is compliant with GHS and OSHA requirements outlined in Hazard Communication Standard 29 CFR 1910.1200. Intended for use only with LDPE and PP compatible solvents.

CATALOG NO.	DESCRIPTION
MLS-370062	UniSafe™ Toluene Vented Wash Bottle
MLS-370061	UniSafe™ Water Vented Wash Bottle
MLS-370060	UniSafe™ Acetonitrile Vented Wash Bottle
MLS-370059	UniSafe™ Sodium Hypochlorite Vented Wash Bottle
MLS-370058	UniSafe™ Methyl Ethyl Ketone Vented Wash Bottle
MLS-370057	UniSafe™ Isopropanol Vented Wash Bottle
MLS-370056	UniSafe™ Ethyl Acetate Vented Wash Bottle
MLS-370055	UniSafe™ Dichloromethane Vented Wash Bottle
MLS-370054	UniSafe™ Distilled Water Vented Wash Bottle
MLS-370053	UniSafe™ Methanol Vented Wash Bottle
MLS-370052	UniSafe™ Ethanol Vented Wash Bottle
MLS-370051	UniSafe™ Acetone Vented Wash Bottle

# BOTTLES

## Reagent Bottles



CATALOG NO.	DESCRIPTION
MLS-33461	Reagent Bottles, Wide Mouth, HDPE, 30 ml, Amber
MLS-33462	Reagent Bottles, Wide Mouth, HDPE, 60 ml, Amber
MLS-33463	Reagent Bottles, Wide Mouth, HDPE, 125 ml, Amber
MLS-33464	Reagent Bottles, Wide Mouth, HDPE, 250 ml, Amber
MLS-33465	Reagent Bottles, Wide Mouth, HDPE, 500 ml, Amber
MLS-33466	Reagent Bottles, Wide Mouth, HDPE, 1000 ml, Amber
<b>These Wide Mouth, HDPE, reagent bottles are manufactured from USP Class VI resin that is FDA-21CFR compliant. Excellent chemical resistance to most acids, bases, and alcohols. Polypropylene caps are included.</b>	

CATALOG NO.	DESCRIPTION
MLS-33421	Reagent Bottles, Narrow Mouth, HDPE, 4 ml, Amber
MLS-33422	Reagent Bottles, Narrow Mouth, HDPE, 8 ml, Amber
MLS-33423	Reagent Bottles, Narrow Mouth, HDPE, 15 ml, Amber
MLS-33424	Reagent Bottles, Narrow Mouth, HDPE, 30 ml, Amber
MLS-33425	Reagent Bottles, Narrow Mouth, HDPE, 60 ml, Amber
MLS-33426	Reagent Bottles, Narrow Mouth, HDPE, 125 ml, Amber
MLS-33427	Reagent Bottles, Narrow Mouth, HDPE, 250 ml, Amber
MLS-33428	Reagent Bottles, Narrow Mouth, HDPE, 500 ml, Amber
MLS-33429	Reagent Bottles, Narrow Mouth, HDPE, 1000 ml, Amber
<b>These narrow Mouth, HDPE, reagent bottles are manufactured from USP Class VI resin that is FDA-21CFR compliant. Excellent chemical resistance to most acids, bases, and alcohols. Polypropylene caps are included.</b>	





# BOTTLES

## Reagent Bottles



CATALOG NO.	DESCRIPTION
MLS-BNMP2L	Reagent Bottles, Narrow Mouth, PP, 2L
MLS-BNMP4L	Reagent Bottles, Narrow Mouth, PP, 4L
MLS-BNMP8L	Reagent Bottles, Narrow Mouth, PP, 8L
<p>These large capacity reagent bottles feature narrow mouths, integrated shoulder loops (except 2L) and include a lined HDPE cap for a tight seal. Constructed of polypropylene, the bottles are translucent, pliable, and have excellent strength for extended use. Ideal for packing, shipping and storage purposes. Cap needs to be removed before autoclaving.</p>	

CATALOG NO.	DESCRIPTION
MLS-BNM090	Reagent Bottles, HDPE, Narrow Mouth, 90 ml
<p>This unique narrow mouth bottle has a capacity of 90ml and a lined polypropylene lid. Constructed of high density polyethylene, the bottle is translucent, pliable, and has excellent strength.</p>	



## Spray Bottles



CATALOG NO.	DESCRIPTION
MLS-SB-8OZ	Spray Bottles with Pump, PET, 8oz
MLS-SB-16OZ	Spray Bottles with Pump, PET, 16oz
<p>Spray bottles are made from clear polyethylene terephthalate (PET) and include a white spray nozzle with clear cap to keep nozzle clean and leak-proof. PET is lightweight, chemically resistant and has excellent durability. Available in 8 oz and 16 oz capacity.</p>	

# BOTTLES

## Dropper Bottle



CATALOG NO.	DESCRIPTION
MLS-PDB015	Dropper Bottle, Unassembled, LDPE, 15ml
<b>This dropper bottle is made of translucent and low density polyethylene. Bottle can be used with or without dropper.</b>	

## Dropping Bottles

CATALOG NO.	DESCRIPTION
MLS-22101	Dropping Bottles Assembled, 60ml
MLS-22102	Dropping Bottles Assembled, 125ml
<b>Dropping bottles dispense small uniform drops of liquids when pressed lightly. Made from low-density polyethylene (LDPE).</b>	



## Centrifuge Bottles



CATALOG NO.	DESCRIPTION
MLS-6330	Centrifuge Bottle, 250mL, Attached Cap, PP
MLS-6332	Centrifuge Bottle, 500mL, Attached Cap, PP
<b>These rugged Polypropylene Centrifuge Bottles are suitable for low to moderate speed centrifugation when filled to 80% of their nominal capacity. Re-usable and autoclavable, the bottles are designed with leakproof linerless closures. Bottles are molded with consistent wall thickness for excellent performance and have excellent chemical resistance.</b>	

# BOTTLES

## Media/Storage Bottles

Clear, square disposable polyethylene terephthalate (PET) bottles are space-saving and convenient to handle. Useful for mixing, sampling, and storage applications. Includes a blue, liner-less polypropylene plug-seal autoclavable cap and features raised graduations. Neck ID is 37mm. PET is transparent, lightweight, and break-resistant. Bottles can be used in environments ranging from -40° to 50° Celsius. Sterilized using gamma radiation. Individually bagged. PET cannot be autoclaved.



CATALOG NO.	DESCRIPTION
MLS-BMPE250	Media/Storage Bottle, Square, PET With PP Cap, 250 ml
MLS-BMPE500	Media/Storage Bottle, Square, PET With PP Cap, 500ML
MLS-BMPE1000	Media/Storage Bottle, Square, PET With PP Cap, 1000ML

Clear, square disposable polycarbonate bottles are space-saving and convenient to handle. Useful for mixing, sampling, and storage applications. Graduations and large marking spot are screened in easy-to-read white enamel. Includes a blue, liner-less GL 45 polypropylene plug-seal autoclavable cap. Neck ID is 37mm. Polycarbonate bottles can be used in environments ranging from -80° to 130° Celsius. Sterilized using gamma radiation. Individually bagged. Bottles can be autoclaved once - repeated autoclaving is not recommended

CATALOG NO.	DESCRIPTION
MLS-BMPC150	Media/Storage Bottle, Square, PC With PP Cap, 150 ml
MLS-BMPC250	Media/Storage Bottle, Square, PC With PP Cap, 250 ml
MLS-BMPC500	Media/Storage Bottle, Square, PC With PP Cap, 500ML
MLS-BMPC1000	Media/Storage Bottle, Square, PC With PP Cap, 1000ML



# ERLENMEYER FLASKS



CATALOG NO.	DESCRIPTION
MLS-FP0125	Wide-mouth Erlenmeyer Flask, PP, 125 ml
MLS-FP0250	Wide-mouth Erlenmeyer Flask, PP, 250 ml
MLS-FP0500	Wide-mouth Erlenmeyer Flask, PP, 500 ml
<b>Wide mouth polypropylene flasks are autoclavable and include a screw cap for easy mixing and storage of solutions. Flasks are graduated to show approximate volumes.</b>	

## VOLUMETRIC FLASK

CATALOG NO.	DESCRIPTION
MLS-BT-V107097	Volumetric Flask, Perfluoxy (PFA) Class A, 10mL
MLS-BT-V107197	Volumetric Flask, Perfluoxy (PFA) Class A, 25mL
MLS-BT-V107297	Volumetric Flask, Perfluoxy (PFA) Class A, 50mL
MLS-BT-V107397	Volumetric Flask, Perfluoxy (PFA) Class A, 100mL
MLS-BT-V107497	Volumetric Flask, Perfluoxy (PFA) Class A, 250mL
<b>These volumetric flasks are highly-transparent, easy-to-clean, and meet Class A tolerances per DIN EN ISO 1042. Flasks are manufactured from PFA (a Perfluoxy copolymer), which is highly resistant to heat and chemical contamination.</b>	



# TISSUE CULTURE FLASK



CATALOG NO.	DESCRIPTION
MLS-TCF-25	50ml, 25cm <sup>2</sup> , Tissue culture flask
MLS-TCF-75	250ml, 75cm <sup>2</sup> , Tissue culture flask
MLS-TCF-175	750ml, 175cm <sup>2</sup> , Tissue culture flask
<b>Tissue Culture Flasks with vent caps, Surface treated with vacuum-gas plasma, Sterilized by gamma radiation</b>	

## FUNNELS

CATALOG NO.	DESCRIPTION
MLS - FPPC65	Funnels, Utility, PP, 65mm
MLS - FPPC52	Funnels, Utility, PP, 52mm
MLS - FPPC35	Funnels, Utility, PP, 35mm
<b>Wide mouth polypropylene flasks are autoclavable and include a screw cap for easy mixing and storage of solutions. Flasks are graduated to show approximate volumes.</b>	





# CARBOYS

## Round Carboys



CATALOG NO.	DESCRIPTION
MLS - 34001	Carboy, PP, 10-Liter
MLS - 34002	Carboy, PP, 20-Liter
<b>These polypropylene, autoclavable carboys feature wide handles for convenient carrying and pouring. Ideal for storing various laboratory solutions. Threads should be completely disengaged before autoclaving.</b>	

## Round Carboys with stopcock

CATALOG NO.	DESCRIPTION
MLS - 34051	Carboy W/Stopcock, PP, 10-Liter
MLS - 34052	Carboy W/Stopcock, PP, 20-Liter
<b>These polypropylene, autoclavable carboys are equipped with a preassembled leak-proof stopcock. Ideal for storing and dispensing various laboratory solutions. Threads should be completely disengaged before autoclaving.</b>	



## Wide Mouth Carboys



CATALOG NO.	DESCRIPTION
MLS - 7260010	Carboys, Round with Handles, Wide Mouth, LDPE, White PP Screwcap, 10 Liter, Molded Graduations
MLS - 7260015	Carboy, Round with Handles, Wide Mouth, LDPE, White PP Screwcap, 15 Liter, Molded Graduations
MLS - 7260020	Carboys, Round with Handles, Wide Mouth, LDPE, White PP Screwcap, 20 Liter, Molded Graduations
<b>Rugged and durable, Carboys are the safe, shatterproof alternative to glass for mixing, storing, and dispensing liquids in volumes from five to fifty liters.</b>	

# CARBOYS

## Heavy Duty Carboys



CATALOG NO.	DESCRIPTION
MLS - 7240010	Carboys, Round with Handles, Heavy Duty PP, White PP Screwcap, 10 Liter, Molded Graduations, Autoclavable
MLS - 7240020	Carboys, Round with Handles, Heavy Duty PP, White PP Screwcap, 20 Liter, Molded Graduations, Autoclavable
<b>These polypropylene, autoclavable carboys feature wide handles for convenient carrying and pouring. Ideal for storing various laboratory solutions. Threads should be completely disengaged before autoclaving.</b>	

## Rectangular Carboys with Spigot

CATALOG NO.	DESCRIPTION
MLS-7300005	Carboys, Rectangular with Spigot and Handle, PP, White PP Screwcap, 5 Liter, Molded Graduations, Autoclavable
MLS-7300010	Carboys, Rectangular with Spigot and Handle, PP, White PP Screwcap, 10 Liter, Molded Graduations, Autoclavable
MLS-7300020	Carboys, Rectangular with Spigot and Handle, PP, White PP Screwcap, 20 Liter, Molded Graduations, Autoclavable
<b>Rugged and durable, Carboys are the safe, shatterproof alternative to glass for mixing, storing, and dispensing liquids in volumes from five to fifty liters.</b>	



## Rectangular Carboys



CATALOG NO.	DESCRIPTION
MLS-601780	Carboys, 2 Liter (1/2 Gallon), HDPE
MLS-601781	Carboys, 5 Liter (1.3 Gallon), HDPE
MLS-601782	Carboys, 10 Liter (2.5 Gallon), HDPE
MLS-601783	Carboys, 20 Liter (5 Gallon), HDPE
<b>Produced from High Density Polyethylene (HDPE), they feature comfortable handles and tamper evident, secure seal caps that follow DIN recommendations.</b>	

# CELL SPREADERS



CATALOG NO.	DESCRIPTION
MLS-HS8171B	10 Spreaders Per Zip Bag; Sterile
MLS-HS8171A	Spreader Paper Peel; Individual Wrap; Sterile
The polypropylene cell spreaders have a smooth surface that will not tear agar or disturb media layers.	

## CELL LIFTERS

CATALOG NO.	DESCRIPTION
MLS-229306	Double Ended Cell Lifter, Flat Blade and J-Hook, Sterile
MLS-229305	Double Ended Cell Lifter, Flat Blade and Narrow Blade, Sterile
Unique double end design features a 19mm beveled edge flat blade on one end, and J-Hook or Narrow Blade on the other	



## CELL SCRAPERS



CATALOG NO.	DESCRIPTION
MLS-229310	Cell Scraper, 25cm Handle, 20mm Scraper Blade, Sterile
MLS-229311	Cell Scraper, 25cm Handle, 20mm Lifter Blade, Sterile
MLS-229312	Cell Scraper, 30cm Handle, 20mm Pivoting Blade, Sterile
MLS-229315	Cell Scraper, 39cm Handle, 30mm Scraper Blade, Sterile
MLS-229316	Cell Scraper, 39cm Handle, 30mm Lifter Blade, Sterile
Cell Scrapers are designed to gently harvest cells from tissue culture flasks, dishes, or bottles	

# CELL STRAINERS



CATALOG NO.	DESCRIPTION
MLS-229480	Pestle for Cell Strainer, Individually Wrapped, Sterile
MLS-229481	40µm Cell Strainer, Individually Wrapped, Sterile
MLS-229482	40µm Cell Strainer, Bulk Packed, Sterile
MLS-229483	70µm Cell Strainer, Individually Wrapped, Sterile
MLS-229484	70µm Cell Strainer, Bulk Packed, Sterile
MLS-229485	100µm Cell Strainer, Individually Wrapped, Sterile
MLS-229486	100µm Cell Strainer, Bulk Packed, Sterile
<p><b>Polypropylene housing with nylon mesh filter eliminated the need to fabricate filters from gauze, paper, nylon mesh, or other materials.</b></p> <p><b>Extended lip on housing enables aseptic handling with forceps</b></p> <p><b>Elongated gripping tab on housing for easy handling</b></p>	

# TEST TUBE PEG RACKS

CATALOG NO.	DESCRIPTION
MLS-456150	Peg Tube Rack, Reinforced PP, 13mm, 96-Place (117 Pegs)
MLS-456155	Peg Tube Rack, Reinforced PP, 17mm, 50-Place (66 Pegs)
<p><b>Fully Autoclavable; Manufactured from durable reinforced polypropylene (PP); Peg design allows easy viewing for tube contents; Can be used to dry inverted tubes while minimizing collection or airborne contaminants.</b></p> <p><b>Available in three colors: Red, White and Blue.</b></p>	



# PCR STRIPS



CATALOG NO.	DESCRIPTION
MLS-PCR-12	PCR 12-well Strips, thin-walled, domed caps
MLS-PCR-08	PCR 8-well Strips, thin-walled, domed caps
MLS-PCRF-08	PCR 8-well Strips, thin-walled, flat caps

# PCR TUBES

CATALOG NO.	DESCRIPTION
MLS-PCR-02	0.2ml PCR Tubes, thin-walled, domed caps
MLS-PCRF-02	0.2ml PCR Tubes, thin-walled, flat caps
MLS-PCRF-05	0.5ml PCR Tubes, thin-walled, flat caps



# WELL PLATES



CATALOG NO.	DESCRIPTION
MLS-WP-06	6-well Plate with lid, individually wrapped, Sterile, flat bottom
MLS-WP-12	12-well Plate with lid, individually wrapped, Sterile, flat bottom
MLS-WP-24	24-well Plate with lid, individually wrapped, Sterile, flat bottom
MLS-WP-96	96-well Plates with lid, individually wrapped, Sterile, Flat bottom
MLS-WPS-96	96-well Plate with 12-well strip, individually wrapped, Sterile, Flat bottom
MLS-WPU-96	96-well Plates with lid, individually wrapped, Sterile, U-shaped bottom
MLS-WPV-96	96-well Plates, no lid, individually wrapped, Sterile, V-shaped bottom
MLS-PCR-96	PCR 96-well Plates, thin-walled, no skirt
MLS-PCRH-96	PCR 96-well Plates, thin-walled, half skirt

# REAGENT RESERVOIRS



CATALOG NO.	DESCRIPTION
MLS-229561	Multiple Cavity Reagent Reservoir, 8 Well, Trough Bottom, PP, Non-sterile
MLS-229562	Multiple Cavity Reagent Reservoir, 12 Well, Trough Bottom, PP, Non-sterile
MLS-229563	Single Cavity Reagent Reservoir, 96 Pyramid Bottom, PP, Non-sterile
MLS-229564	Single Cavity Reagent Reservoir, 384 Pyramid Bottom, PP, Non-sterile
MLS-229565	Multiple Cavity Reagent Reservoir, 8 Well, Trough Bottom, PP, Sterile
MLS-3054-1010	Multiple Cavity Reagent Reservoir, 12 Well, Trough Bottom, PP, Sterile
MLS-3054-1011	12 Channel Reagent Reservoir, Individually Wrapped, Sterile
MLS-4054-1010	12 Channel Reagent Reservoir, Individually Wrapped, Sterile
MLS-4054-1011	12 Channel Reagent Reservoir, Bulk, Non-sterile
MLS-4054-1011	12 Channel Reagent Reservoir, Bulk, Non-sterile
<p><b>Reservoirs are designed to assist in repetitive or robotic liquid transfer. Pyramid and Trough Bottom designs promote funneling of liquid to maximize reagent recovery by pipette tips. Low profile is ideal for small volume robotic tips. Can be stored at temperatures as low as -80°C. Polypropylene offers a low binding surface to prevent samples from sticking to the side walls during elution, and is chemically inert for combinatorial chemistry applications.</b></p>	



# POLYCARBONATE BIOBOX



CATALOG NO.	DESCRIPTION
MLS-3039	BioBox 25, for 1.0mL and 2.0mL CryoClear Vials, Polycarbonate (PC), Holds 25 vials (5x5 format), Printed Lid, Pack Includes a CryoClear Vial Picker
Available in four colors: Blue, Green, Red, Yellow.	

## FREEZER STORAGE BOXES

CATALOG NO.	DESCRIPTION
MLS-3057	Freezer Box for 5.0mL Microcentrifuge Tubes, 21-place (7 x 3 format), with Lid, Clear
MLS-3055	Freezer Box (no wells or slots), 139.7 x 71.1 x 53.0mm
MLS-3056	Freezer Box for 1.5mL and 2.0mL Microcentrifuge Tubes, 50-place (5 x 10 format), with Lid, Clear
These rugged polypropylene boxes are designed for the compact storage of microcentrifuge tubes, 5ml midi tubes, sample tubes, and cryogenic vials in a temperature range of -196° to +121°C.	



## ICE BUCKETS WITH LID



CATALOG NO.	DESCRIPTION
MLS-455015	Ice Bucket with Cover, 4.5 Liter
Rectangular ice buckets are suitable for use with liquid nitrogen, dry ice, wet ice, and salt solutions. Lids included with ice buckets and select ice trays limit evaporation for extended chilling. Molded handles provide a secure grip for easy transport. Available in four colors: Blue, Red, Green and Black.	

# WIRELESS NYLON RACKS

These fully autoclavable racks are molded from reinforced nylon for strength and durability. With no coatings to peel or chip, they are a great alternative to wire coated racks where the formation of rust can be an issue. Seven colorful options to keep samples and benchtop organized.

## Full Size Wireless Racks



CATALOG NO.	DIAMETER	PLACES	ARRAY	DIMENSIONS
MLS-456440	13mm	72	6x12	200 x 102 x 57mm
MLS-456540	16mm	72	6x12	248 x 127 x 70mm
MLS-456640	20mm	40	4x10	250 x 102 x 83mm
MLS-456740	25mm	40	4x10	300 x 121 x 92mm
MLS-456840	30mm	24	3x8	283 x 108 x 83mm

## Half Size Wireless Racks

CATALOG NO.	DIAMETER	PLACES	ARRAY	DIMENSIONS
MLS-456445	13mm	36	6 x 6	102 x 102 x 56mm
MLS-456545	16mm	36	6 x 6	126 x 126 x 68mm
MLS-456645	20mm	20	4 x 5	128 x 103 x 83mm
MLS-456745	25mm	16	4 x 4	122 x 122 x 75mm
MLS-456845	30mm	9	3 x 3	109 x 109 x 84mm



## Wireless Racks for Microtubes



CATALOG NO.	DESCRIPTION
MLS-456345	Wireless MicroTube Rack with Handles for 1.5mL and 2.0mL Microcentrifuge Tubes, 24-Place
MLS-456340	Wireless MicroTube Rack with Handles for 1.5mL and 2.0mL Microcentrifuge Tubes, 96-Place
Four colorful options to keep samples and benchtop organized	

# STIR BARS

## Cylindrical Stir Bars



CATALOG NO.	DESCRIPTION
MLS-SWN652	Stir Bars, Cylindrical, PTFE, 6mm
MLS-SWN654	Stir Bars, Cylindrical, PTFE, 8mm
MLS-SWN656	Stir Bars, Cylindrical, PTFE, 12mm
MLS-SWN658	Stir Bars, Cylindrical, PTFE, 15mm
MLS-SWN660	Stir Bars, Cylindrical, PTFE, 20mm
MLS-SWN662	Stir Bars, Cylindrical, PTFE, 25mm
<b>Strong alloy magnets sheathed in pure white inert PTFE. Resistant to virtually all chemical attack and temperatures up to 260°C .</b>	

CATALOG NO.	DESCRIPTION
MLS-SWN664	Stir Bars, Cylindrical, PTFE, 30mm
MLS-SWN666	Stir Bars, Cylindrical, PTFE, 35mm
MLS-SWN668	Stir Bars, Cylindrical, PTFE, 40mm
MLS-SWN670	Stir Bars, Cylindrical, PTFE, 45mm
MLS-SWN672	Stir Bars, Cylindrical, PTFE, 50mm
MLS-SWN674	Stir Bars, Cylindrical, PTFE, 60mm
<b>Strong alloy magnets sheathed in pure white inert PTFE. Resistant to virtually all chemical attack and temperatures up to 260°C .</b>	

## Micro Stir Bars



CATALOG NO.	DESCRIPTION
MLS-SWN950	Stir Bars, Micro, PTFE, 7mm
MLS-SWN954	Stir Bars, Micro, PTFE, 10mm
MLS-SWN956	Stir Bars, Micro, PTFE, 13mm
<b>Micro or flea magnetic stir bars are designed for stirring small volumes in vessels such as vials, tubes and gradient makers.</b>	

# STIR BARS

## Oval Stir Bars



CATALOG NO.	DESCRIPTION
MLS-SWN520	Stir Bars, Oval, PTFE, 20mm
MLS-SWN522	Stir Bars, Oval, PTFE, 25mm
MLS-SWN524	Stir Bars, Oval, PTFE, 30mm
MLS-SWN525	Stir Bars, Oval, PTFE, 35mm
MLS-SWN526	Stir Bars, Oval, PTFE, 40mm
MLS-SWN528	Stir Bars, Oval, PTFE, 50mm
<b>Produced from PTFE For efficient stirring in round bottomed flasks</b>	

## Pivot Ring Stir Bars



CATALOG NO.	DESCRIPTION
MLS-SWN501	Stir Bars, Pivot Ring, PTFE, 15mm
MLS-SWN502	Stir Bars, Pivot Ring, PTFE, 12mm
MLS-SWN503	Stir Bars, Pivot Ring, PTFE, 20mm
MLS-SWN504	Stir Bars, Pivot Ring, PTFE, 25mm
MLS-SWN505	Stir Bars, Pivot Ring, PTFE, 35mm
<b>Produced from PTFE. Central pivot ring ensures efficient spinning, even on curved or uneven bases.</b>	

CATALOG NO.	DESCRIPTION
MLS-SWN506	Stir Bars, Pivot Ring, PTFE, 45mm
MLS-SWN507	Stir Bars, Pivot Ring, PTFE, 40mm
MLS-SWN508	Stir Bars, Pivot Ring, PTFE, 60mm
MLS-SWN509	Stir Bars, Pivot Ring, PTFE, 50mm
MLS-SWN510	Stir Bars, Pivot Ring, PTFE, 70mm
<b>Produced from PTFE. Central pivot ring ensures efficient spinning, even on curved or uneven bases.</b>	

# WEIGHING BOATS

## Hexagonal Weighing Boats



CATALOG NO.	DESCRIPTION
MLS-3615	Weighing Dish, Plastic, Hexagonal, Antistatic, 20mL, PS
MLS-3616	Weighing Dish, Plastic, Hexagonal, Antistatic, 50mL, PS
MLS-3617	Weighing Dish, Plastic, Hexagonal, Antistatic, 200mL, PS
MLS-3618	Weighing Dish, Plastic, Hexagonal, Antistatic, 350mL, PS
Hexagonal Weighing Boats are uniquely designed for non-slip handling, pouring, and dispensing of solids and liquids when weighing. The weigh boats have flat bottoms and sloping sides to ensure stability, prevent tipping and facilitate accurate pour-outs.	

## Diamond Shaped Weighing Boats



CATALOG NO.	DESCRIPTION
MLS-3614-100	Diamond Shaped Weighing Boat, Anti-static, 125 x 100 x 18mm, 100ml
MLS-3614-30	Diamond Shaped Weighing Boat, Anti-static, 80 x 60 x 18mm, 30ml
MLS-3614-5	Diamond Shaped Weighing Boat, Anti-static, 58 x 36 x 6mm, 5ml
Diamond Shaped Weighing Boats are designed for a wide range of applications including weighing, dispensing, and storing of liquids and solids. The flexible weigh boats have a smooth surface and sloping sides that provide controlled and accurate pour-outs with minimal sample loss.	

# WEIGHING BOATS

## Square Weighing Boats



CATALOG NO.	DESCRIPTION
MLS-3619-10	Weight Boat, Square with Square Bottom, Antistatic, PS, White, 10mL
MLS-3619-100	Weight Boat, Square with Square Bottom, Antistatic, PS, White, 100mL
MLS-3619-250	Weight Boat, Square with Square Bottom, Antistatic, PS, White, 250mL
MLS-3619-10BK	Weight Boat, Square with Square Bottom, Antistatic, PS, Black, 10mL
MLS-3619-100BK	Weight Boat, Square with Square Bottom, Antistatic, PS, Black, 100mL
MLS-3619-250BK	Weight Boat, Square with Square Bottom, Antistatic, PS, Black, 250mL
MLS-3620	Weighing Boat, Plastic, Square with Round Bottom, Small Easy Pour Spout, Antistatic, 41 x 41 x 8mm, PS, White, 7mL
MLS-3621	Weighing Boat, Plastic, Square with Round Bottom, Small Easy Pour Spout, Antistatic, 89 x 89 x 25mm, PS, White, 100mL
MLS-3622	Weighing Boat, Plastic, Square with Round Bottom, Small Easy Pour Spout, Antistatic, 140 x 140 x 25mm, PS, White, 330mL
<b>Weighing Boats are designed for a wide range of applications including weighing, dispensing, and storing of liquids and solids. They can also be used as quick freeze trays for sample material, discard trays for broken ampoules, or mixing trays for small batches. The flexible weigh boats have a smooth surface, sloping sides and recessed rounded corners that provide controlled and accurate pour-outs with minimal sample loss.</b>	



# ALUMINUM DISHES

## Round Aluminum Dishes with Tabs



CATALOG NO.	DIAMETER	CAPACITY
MLS-D28-500	28mm	8ml
MLS-D43-100	43mm	20ml
MLS-D57-100	57mm	50ml
MLS-D70-100	70mm	60ml

## Smooth-Walled Aluminum Weighing Dishes

CATALOG NO.	DIAMETER	CAPACITY
MLS-D57S-100	57mm	50ml
MLS-D70S-100	73mm	80ml
MLS-D140S-100	103mm	140ml



## General Purpose Aluminum Weighing Dishes



CATALOG NO.	DIAMETER	CAPACITY
MLS-D44-100	44mm	20ml
MLS-D63-100	63mm	70ml
MLS-D75-100	62mm	75ml
MLS-D200-50	110mm	200ml
MLS-D500-50	181mm	480ml
MLS-111	95mm	250ml

# ALUMINUM DISHES

## Aluminum Weighing/Drying Pans



CATALOG NO.	DIAMETER	WEIGHT
MLS-D123-50	70mm	1.2 g
MLS-D126-50	102mm	2.5 g
MLS-D125-100	127mm	3.3 g
MLS-D11-25	110mm	3.7 g
MLS-D12-25	110mm	4.4 g

## Micro Aluminum Weighing Dishes



CATALOG NO.	SIZE	CAPACITY
MLS-MAWD-0400	6.5mm x 3.5mm x 2.5mm	0.4 ml
MLS-MAWD-1100	10mm x 4mm x 4.5mm	0.11 ml
MLS-MAWD-1200	12mm x 4mm x 2.7mm	0.12 ml
MLS-MAWD-1201	12mm x 2mm x 4.5mm	0.12 ml
MLS-MAWD-1120	25mm x 7mm x 7mm	1.12 ml
MLS-MAWD-3500	13mm x 3.5mm	0.35 ml
MLS-MAWD-0350	6mm x 2.5mm	0.35 ml

# PLASTIC TRAYS

## Pipette/Instrument Organizer Tray



CATALOG NO.	DESCRIPTION
MLS-PPT004	Pipette/Instrument Organizer Tray, PP
Tray features 4 deep columns (approx.. ¾" deep) and is recommended for 1 to 10ml pipettes. Dimensions are 16" x 12" x 1" and is autoclavable. Molded in polypropylene and features 1/8" thick walls.	

## Laboratory Trays

CATALOG NO.	DESCRIPTION
MLS-81722	Laboratory tray, Large, PP, 20" X 17" X 5"
MLS-81701	Laboratory tray, PP, 18" X 14" X 3"
MLS-81702	Laboratory tray, PP, 18" X 14" X 3"
MLS-81721	Utility Trays, PP, 15" X 14" X 5"
Large polypropylene molded laboratory tray is steam autoclavable and can be used for sterilizing and drying labware.	



## Draining Tray



CATALOG NO.	DESCRIPTION
MLS-78103	Draining Tray, PP, 16" X 16" X 4"
This polypropylene basket with wide net-like structure on all sides is useful for post-wash draining of laboratory apparatus. Autoclavable	

## Carrier Tray

CATALOG NO.	DESCRIPTION
MLS-81731	Carrier Tray, PP, 15" X 9.5" X 4.5"
Polypropylene molded carrier tray is steam autoclavable. Includes a builtin handle for easy carrying of reagent bottles or other laboratory ware using one hand.	



# STANDS

## Micropipette Stand



CATALOG NO.	DESCRIPTION
MLS-MPSM01	Micropipette Stand, PP
<b>Ideal for conveniently holding micropipettes at a workstation. This polypropylene stand is designed for easy access and features a storage compartment for tips.</b>	

## Flask Stand

CATALOG NO.	DESCRIPTION
MLS-45105	Flask Stand
<b>For use with flasks up to 10L, this polypropylene stand is ribbed to safely hold flasks in place. Autoclavable.</b>	



## Pipette Stand



CATALOG NO.	DESCRIPTION
MLS-79101	Pipette Stand, Horizontal, 12-Place, PP
<b>Polypropylene pipette stand can hold up to 12 pipettes in a horizontal position. Autoclavable.</b>	

## Pipette Stand

CATALOG NO.	DESCRIPTION
MLS-79102	Pipette Stand, Vertical, 28-Place, PP
<b>Polypropylene, non-corrosive pipette stand can hold up to 28 pipettes. The heavy base is 9" in diameter and the polypropylene coated square rod is 15" high.</b>	



CATALOG NO.	DESCRIPTION
MLS-79103	Pipette Stand, Rotary, 94-Place, PP
<b>Central vertical axis measures 19" in height. The lower rotating disc measures 9" in diameter and has four concentric troughs that align with the upper disc.</b>	



# CUVETTES



CATALOG NO.	DESCRIPTION
MLS-CUV-04	Macro Four sided Cuvettes, Clear
<p><b>cuvettes are manufactured from high quality polystyrene and formulated for superior clarity and transmission. These fit most spectrophotometers including Beckman, Perkin-Elmer, Gilford, Hitachi, Varian, Shimadzu and are packed in dust resistant trays.</b></p>	

CATALOG NO.	DESCRIPTION
MLS-CUV-112117	Cuvette, Semi-Micro, 2.9mL, with 2 Clear Sides, PS
<p><b>These cuvettes are engineered for high reproducibility and offer variations within the range of <math>\pm 1\%</math> absorption. This is vital when batch analysis is being undertaken.</b></p>	



CATALOG NO.	DESCRIPTION
MLS-CUV-112137	Cuvette, Semi-Micro, 2.9mL, with 2 Clear Sides, PS
<p><b>These cuvettes are engineered for high reproducibility and offer variations within the range of <math>\pm 1\%</math> absorption. This is vital when batch analysis is being undertaken.</b></p>	

CATALOG NO.	DESCRIPTION
MLS-CUV-111137	Cuvette, Spectrophotometer, Square, 4.5mL (10mm), PS, 4 Clear Sides
<p><b>These cuvettes are engineered for high reproducibility and offer variations within the range of <math>\pm 1\%</math> absorption. This is vital when batch analysis is being undertaken.</b></p>	



CATALOG NO.	DESCRIPTION
MLS-CUV-111117	Cuvette, Spectrophotometer, Square, 4.5mL (10mm), PS, 2 Clear Sides
<p><b>These cuvettes are engineered for high reproducibility and offer variations within the range of <math>\pm 1\%</math> absorption. This is vital when batch analysis is being undertaken.</b></p>	

# TRANSPORT TUBE



CATALOG NO.	DESCRIPTION
MLS-111F4	5ml Vol, Multi-Vial, Without Cap
<b>Multi-vials are an excellent multipurpose container, and are primarily used for liquid scintillation counting, beta-gamma counting, and sample storage. These vials include white polyethylene caps or available without caps.</b>	

CATALOG NO.	DESCRIPTION
MLS-55F2	12ml Transport Tube, Polypropylene Material, No Cap, Flat Bottom
MLS-53F2	12ml Transport Tube, Polypropylene Material, No Cap, Conical Skirted Bottom
<b>These 12ml Transport Tubes are made from high quality polypropylene or polystyrene. Our Transport Tubes are ideal for transporting and collecting specimens. These tubes have a flat bottom.</b>	



# URINE CENTRIFUGE TUBES



CATALOG NO.	DESCRIPTION
MLS-112011	Tube, Urine Centrifuge, 12mL, with Flared Top, PS, Graduated
MLS-112014	Urine Collection System, 12mL Flared Top Urine Centrifuge Tube and Separate Blue Snap Cap
MLS-112017	Tube, Urine Centrifuge, 12mL, with Flared Top, PS, Graduated
MLS-112030	Tube, Urine Centrifuge, 12mL, with Sediment Bulb (Whale Type) and Flared Top, PS, Graduated to 10 mL
MLS-112015	Tube, Urine Centrifuge, 12mL with Flared Top, PS, Kova Type, Graduated
<b>This line of centrifuge tubes has been designed specifically for urinalysis procedures. The tubes offer a flared-opening which is designed to make pouring off easy.</b>	



# TEST TUBES



CATALOG NO.	DESCRIPTION
MLS-110410	Test Tube, 12 x 75mm, (5mL), PS
MLS-110441	Test Tube, 12 x 75mm, (5mL), PP
MLS-110416	Test Tube, 12 x 75mm, (5mL), PS
MLS-110446	Test Tube, 12 x 75mm, (5mL), PP
MLS-110409	Test Tube, 12 x 75mm, (5mL), PS
MLS-110440	Test Tube, 12 x 75mm, (5mL), PP
MLS-110410LB	Test Tube, 12 x 75mm, (5mL), PS, Light Blue
MLS-110441AM	Test Tube, 12 x 75mm, (5mL), PP, Amber
MLS-110413	Test Tube, 12 x 75mm, (5mL), PS
MLS-110414	Test Tube, 12 x 75mm, (5mL), PS
MLS-110443	Test Tube, 12 x 75mm, (5mL), PP
MLS-110444	Test Tube, 12 x 75mm, (5mL), PP
<b>Plastic culture tubes are produced in polystyrene (PS) or polypropylene (PP). The polystyrene tubes can withstand 1500 RPM during centrifugation. polypropylene tubes can withstand 3000 RPM and are autoclavable and freezable.</b>	

# PITCHERS



CATALOG NO.	DESCRIPTION
MLS-601156	Pitcher, PP, Molded Graduations, 500mL
MLS-601157	Pitcher, PP, Molded Graduations, 1000mL
MLS-601158	Pitcher, PP, Molded Graduations, 2000mL
MLS-601159	Pitcher, PP, Molded Graduations, 3000mL
MLS-601160	Pitcher, PP, Molded Graduations, 5000mL
<b>They are made from heavy-duty plastic and withstands repeated autoclaving.</b>	

# CRYO VIALS

CATALOG NO.	DESCRIPTION
MLS-P60109	Cryo vial, Internal Thread, 1.8ml, Star Base
MLS-P60110	Cryo vial, Internal Thread, 1.8ml Round Base
MLS-P60112	Cryo vial, Internal Thread, 4.5ml, Star Base
MLS-P60113	Cryo vial, Internal Thread, 4.5ml Round Base
<b>Cryo vials are designed for storage of samples as low as -196°C*. DNase, RNase, and Pyrogen-free.</b>	



# CRYO STORAGE BOX



CATALOG NO.	DESCRIPTION
MLS-P20607	Hinged Storage Box for Microtubes/Cryo vials, 100 places
<b>Cryo storage boxes are designed for storage of samples as low as -90°C. Hold both internal and external threaded cryogenic vials as well as storage vials.</b>	

# SPATULAS AND SCOOPS

## Metric Measuring Spoon Set



CATALOG NO.	DESCRIPTION
MLS-UNSPOON-B	Metric Measuring Spoon Set
<b>Unbreakable plastic set of five metric measuring spoons. Includes 1.25ml (1/4 tsp), 2.5ml (1/2 tsp), 5ml (1 tsp), 7.5ml (1/2 tb sp) and 15ml (1 tb sp).</b>	

## Laboratory Scoops

CATALOG NO.	DESCRIPTION
MLS-81253	Laboratory Scoops, PP, 10 ml
MLS-81254	Laboratory Scoops, PP, 25 ml
MLS-81255	Laboratory Scoops, PP, 50 ml
MLS-81256	Laboratory Scoops, PP, 100 ml
MLS-81257	Laboratory Scoops, PP, 250 ml
MLS-81258	Laboratory Scoops, PP, 500 ml
MLS-81259	Laboratory Scoops, PP, 100 ml
<b>These sturdy polypropylene scoops are available in a variety of sizes and as a complete set of seven scoops. Capacities are molded into the top of each scoop for easy reference. The PP material is autoclavable and will not rust or corrode.</b>	




## Spatulas



CATALOG NO.	DESCRIPTION
MLS-81411	Spatula, 9.5" long, 1/4" Diameter
MLS-81301	spatulas feature a spoon on one side, 6"
MLS-81302	spatulas feature a spoon on one side, 8"



# PLASTIC TWEEZERS



CATALOG NO.	DESCRIPTION
MLS-PI012-1	Plastic Tweezers
<b>These polypropylene tweezers are 4" long and feature serrated tips. Non-sterile, autoclavable</b>	

# PLASTIC FORCEPS

CATALOG NO.	DESCRIPTION
MLS-FOBPP5	Plastic Forceps, 4.5"
<b>These 4.5" long ABS plastic forceps feature serrated inside tips and ribbed finger supports. Non-sterile</b>	




# TISSUE CASSETTES



CATALOG NO.	DESCRIPTION
MLS-1094	Cassette, Tissue Embedding with Attached Lid, 35° Writing Area
<b>The lid securely locks into the base and keeps samples safe and secure. The 35° angled cassettes can be manually loaded into most cassette printers. Available in ten different shades.</b>	

# CAPILLARY TUBES

CATALOG NO.	DESCRIPTION
MLS-51680	Capillary Tube, Microhematocrit, 100% Plastic, Blue Tip, Plain
MLS-51684	Capillary Tube, Microhematocrit, 100% Plastic, Green Tip, Ammonium Heparinized
MLS-51686	Capillary Tube, Microhematocrit, 100% Plastic, Red Tip, Sodium Heparinized



# PACK SIZES

## Centrifuge Tubes

CATALOG NO.	PACK SIZE
MCT-00015	50 pcs/pack
MCT-00050	50 pcs/pack
MLS-3131-345-008	50 pcs/pack
MLS-3181-345-008	50 pcs/pack
MLS-D1002	50 pcs/pack
MLS-3095-340-008	50 pcs/pack
MLS-3144-335-001	50 pcs/pack
MLS-3184-345-001	50 pcs/pack
MLS-10482	500 pcs/pack
MLS-10485	500 pcs/pack
MLS-10480	500 pcs/pack
MLS-10483	500 pcs/pack

## Plastic Funnels

CATALOG NO.	PACK SIZE
MLS - FPPC65	12 pcs/pack
MLS - FPPC52	12 pcs/pack
MLS - FPPC35	12 pcs/pack

## Petri Dishes

CATALOG NO.	PACK SIZE
MPD-00060-500	10 pc/ pack
MPD-00100-500	40 pc/ pack
MPD-00150-500	10 pc/ pack

## Sample Container

CATALOG NO.	PACK SIZE
MLS-SC-120	200 pc/ pack
MLS-C-11008-N	500 pc/ pack

## Serological Pipets

CATALOG NO.	PACK SIZE
MSP-00001	500 pcs/pack
MSP-00002	400 pcs/pack
MSP-00005	200 pcs/pack
MSP-00010	200 pcs/pack
MSP-00025	200 pcs/pack
MSP-00050	25 pcs/pack

## Pipette Filter Tips

CATALOG NO.	PACK SIZE
MLS-FPT-0010	96/rack
MLS-FPT-0020	96/rack
MLS-FPT-0200	96/rack
MLS-FPT-1000	96/rack
MLS-FPT-1250	96/rack

## Incoculating Loops

CATALOG NO.	PACK SIZE
MPD-00060-500	10 pc/ pack
MPD-00100-500	40 pc/ pack

## Racked Pipette Tips

CATALOG NO.	PACK SIZE
MLS-RPT-0010	96/rack
MLS-RPT-0200	96/rack
MLS-RPT-1000	96/rack
MLS-RPT-1250	96/rack

## Cell Spreaders

CATALOG NO.	PACK SIZE
MLS-HS8171B	500/ pack
MLS-HS8171A	500/ pack

## Cell Lifters

CATALOG NO.	PACK SIZE
MLS-229306	1000/ pack
MLS-229305	1000/ pack

## Pipette Tips

CATALOG NO.	PACK SIZE
MPT-00010N	1000 pcs/pack
MPT-00200N	1000 pcs/pack
MPT-01000N	500 pcs/pack
MLS-PT1250	1000 pcs/pack

Looking for a different option? We can do it!  
Contact us now for more details.

# PACK SIZES

## Micro Centrifuge Tubes

CATALOG NO.	PACK SIZE
MLS-MCT1.5	500 pcs/pack
MLS-111642	1000 pcs/pack
MLS-111620	1000 pcs/pack
MLS-111550	1000 pcs/pack
MLS-110550	1000 pcs/pack
MLS-111568	1000 pcs/pack
MLS-111552LK	500 pcs/pack
MLS-111562LK	500 pcs/pack
MLS-111572LK	500 pcs/pack
MLS-111582LK	100 pcs/pack
MLS-110540	1000 pcs/pack
MLS-111554A	500 pcs/pack
MLS-111564A	500 pcs/pack

CATALOG NO.	PACK SIZE
MLS-111574A	500 pcs/pack
MLS-111554B	500 pcs/pack
MLS-111554G	500 pcs/pack
MLS-111554R	500 pcs/pack
MLS-111554Y	500 pcs/pack
MLS-111564B	500 pcs/pack
MLS-111564G	500 pcs/pack
MLS-111564R	500 pcs/pack
MLS-111564Y	500 pcs/pack
MLS-111574B	500 pcs/pack
MLS-111574G	500 pcs/pack
MLS-111574R	500 pcs/pack
MLS-111574Y	500 pcs/pack

## Beakers

CATALOG NO.	PACK SIZE
MLS-PMP-50	12 pcs/pack
MLS-PMP-100	12 pcs/pack
MLS-PMP-250	12 pcs/pack
MLS-PMP-500	6 pcs/pack
MLS-PMP-1000	4 pcs/pack
MLS-PMP-2000	4 pcs/pack
MLS-PMP-5000	4 pcs/pack
MLS-PMP-10000	1 pc/pack

## Beakers with Handle

CATALOG NO.	PACK SIZE
MLS-PBH-100	6 pcs/pack
MLS-PBH-250	6 pcs/pack
MLS-PBH-500	6 pcs/pack
MLS-PBH-1000	4 pcs/pack
MLS-PBH-2000	4 pcs/pack
MLS-PBH-5000	2 pcs/pack
MLS-PBH-10000	1 pc/pack

## Wide Spout Beakers

CATALOG NO.	PACK SIZE
MLS-PP-0025	12 pcs/pack
MLS-PP-0050	12 pcs/pack
MLS-PP-0100	12 pcs/pack
MLS-PP-0250	6 pcs/pack

CATALOG NO.	PACK SIZE
MLS-PP-0500	4 pcs/pack
MLS-PP-1000	3 pcs/pack
MLS-PP-2000	3 pcs/pack
MLS-PP-5000	1 pc/pack

## Transport Tube

CATALOG NO.	PACK SIZE
MLS-111F4	1000 pcs/pack
MLS-55F2	500 pcs/pack
MLS-53F2	500 pcs/pack

Looking for a different option? We can do it!  
Contact us now for more details.



# PACK SIZES

## Griffin Style Beakers

CATALOG NO.	PACK SIZE
MLS-PPG-0050	12 pcs/pack
MLS-PPG-0100	12 pcs/pack
MLS-PPG-0250	6 pcs/pack
MLS-PPG-0500	4 pcs/pack
MLS-PPG-0600	4 pcs/pack
MLS-PPG-1000	3 pcs/pack
MLS-PPG-2000	3 pcs/pack

## Tri-Corner Beakers

CATALOG NO.	PACK SIZE
MLS-PPT-0050	100 pcs/pack
MLS-PPT-0100	100 pcs/pack
MLS-PPT-0250	100 pcs/pack
MLS-PPT-0400	100 pcs/pack
MLS-PPT-0800	100 pcs/pack
MLS-PPT-1000	100 pcs/pack

CATALOG NO.	PACK SIZE
MLS-CPP-0050	100 pcs/pack
MLS-CPP-0100	100 pcs/pack
MLS-CPP-0250	100 pcs/pack
MLS-CPP-0400	100 pcs/pack
MLS-CPP-0800	100 pcs/pack
MLS-CPP-1000	100 pcs/pack

## Erlenmeyer Flasks

CATALOG NO.	PACK SIZE
MLS-FP0125	12 pcs/pack
MLS-FP0250	12 pcs/pack
MLS-FP0500	12 pcs/pack

## Tissue culture flask

CATALOG NO.	PACK SIZE
MLS-TCF-25	200 pcs/pack
MLS-TCF-75	200 pcs/pack
MLS-TCF-175	200 pcs/pack

## cuvettes

CATALOG NO.	PACK SIZE
MLS-CUV-04	100 pcs/pack
MLS-CUV-112117	500 pcs/pack
MLS-CUV-112137	500 pcs/pack
MLS-CUV-111137	500 pcs/pack
MLS-CUV-111117	500 pcs/pack

## Cell Strainers

CATALOG NO.	PACK SIZE
MLS-229480	1 pcs/pack
MLS-229481	1 pc/pack
MLS-229482	10 pcs/pack
MLS-229483	1 pc/pack
MLS-229484	10 pcs/pack
MLS-229485	1 pcs/pack
MLS-229486	10 pcs/pack

## Volumetric Flask

CATALOG NO.	PACK SIZE
MLS-BT-V107097	6 pcs/pack
MLS-BT-V107197	6 pcs/pack
MLS-BT-V107297	6 pcs/pack
MLS-BT-V107397	6 pcs/pack
MLS-BT-V107497	6 pcs/pack

## Urine Centrifuge Tubes

CATALOG NO.	PACK SIZE
MLS-112011	500 pcs/pack
MLS-112014	1000 pcs/pack
MLS-112017	1500 pcs/pack
MLS-112030	1500 pcs/pack
MLS-112015	2000 pcs/pack

**Looking for a different option? We can do it!**  
**Contact us now for more details.**

# PACK SIZES

## Measuring Cylinders

CATALOG NO.	PACK SIZE
MLS-CPSET7-H	Set of 7
MLS-CPSET5	set of 5
MLS-81150	12 pcs/pack
MLS-81151	12 pcs/pack
MLS-81152	12 pcs/pack
MLS-81153	12 pcs/pack
MLS-81154	12 pcs/pack
MLS-P50501	1 pc/pack
MLS-P50502	1 pc/pack
MLS-P50503	1 pc/pack
MLS-P50504	1 pc/pack
MLS-P50505	1 pc/pack
MLS-P50506	1 pc/pack
MLS-P50507	1 pc/pack
MLS-P50508	1 pc/pack

## Pitchers

CATALOG NO.	PACK SIZE
MLS-601156	60 pcs/case
MLS-601157	48 pcs/pack
MLS-601158	36 pcs/pack
MLS-601159	27 pcs/pack
MLS-601160	12 pcs/pack

CATALOG NO.	PACK SIZE
MLS-P50401	12 pcs/pack
MLS-P50402	12 pcs/pack
MLS-P50403	12 pcs/pack
MLS-P50404	6 pcs/pack
MLS-P50405	6 pcs/pack
MLS-P50406	3 pcs/pack
MLS-P50407	3 pcs/pack
MLS-P50408	3 pcs/pack
MLS-80031	24 pcs/pack
MLS-80032	12 pcs/pack
MLS-80033	12 pcs/pack
MLS-80034	12 pcs/pack
MLS-80035	6 pcs/pack
MLS-80036	4 pcs/pack
MLS-80037	3 pcs/pack
MLS-80038	2 pcs/pack

## Test Tubes

CATALOG NO.	PACK SIZE
MLS-110410	1000 pcs/pack
MLS-110441	1000 pcs/pack
MLS-110416	1000 pcs/pack
MLS-110446	1000 pcs/pack
MLS-110409	2000 pcs/pack
MLS-110440	2000 pcs/pack
MLS-110410LB	1000 pcs/pack
MLS-110441AM	2000 pcs/pack
MLS-110413	1000 pcs/pack
MLS-110414	2000 pcs/pack
MLS-110443	1000 pcs/pack
MLS-110444	2000 pcs/pack

## Cryo Vials

CATALOG NO.	PACK SIZE
MLS-P60109	500 pcs/pack
MLS-P60110	500 pcs/pack
MLS-P60112	250 pcs/pack
MLS-P60113	250 pcs/pack

**Looking for a different option? We can do it!  
Contact us now for more details.**

# PACK SIZES

## Transfer Pipets

CATALOG NO.	PACK SIZE
MLS-134010	4000 pcs/pack
MLS-134020	3200 pcs/pack
MLS-134070	2000 pcs/pack
MLS-134050	5000 pcs/pack
MLS-134060	5000 pcs/pack
MLS-134090	4000 pcs/pack
MLS-P31205	500 pcs/pack
MLS-138020	5000 pcs/pack
MLS-138060	5000 pcs/pack
MLS-138050	5000 pcs/pack
MLS-138080	5000 pcs/pack

CATALOG NO.	PACK SIZE
MLS-138030	4000 pcs/pack
MLS-138090	2500 pcs/pack
MLS-138040	5000 pcs/pack
MLS-136030	5000 pcs/pack
MLS-139020	5000 pcs/pack
MLS-139040	2500 pcs/pack
MLS-136036	5000 pcs/pack
MLS-137035	5000 pcs/pack
MLS-138070	5000 pcs/pack
MLS-137050	5000 pcs/pack
MLS-135030	5000 pcs/pack

## Nylon Syringe Filters

CATALOG NO.	PACK SIZE
MSF-NY2213	100 pcs/pack
MSF-NY4513	100 pcs/pack
MSF-NY2225	100 pcs/pack
MSF-NY4525	100 pcs/pack

## Hydrophobic PTFE Syringe Filters

CATALOG NO.	PACK SIZE
MSF-PT2213	100 pcs/pack
MSF-PT4513	100 pcs/pack
MSF-PT2225	100 pcs/pack
MSF-PT4525	100 pcs/pack

## Vent Filters With PTFE Membranes

CATALOG NO.	PACK SIZE
MSF-VF-02713	10 pcs/pack
MSF-VF-02612	10 pcs/pack
MSF-VF-04550	10 pcs/pack

## PES Syringe Filters

CATALOG NO.	PACK SIZE
MSF-PES-0825	500 pcs/pack
MSF-PES-4528	50 pcs/pack
MSF-PES-2028	50 pcs/pack
MSF-PES-2213	100 pcs/pack
MSF-PES-4513	50 pcs/pack
MSF-PES-2217	100 pcs/pack

## Hydrophobic PVDF Syringe Filters

CATALOG NO.	PACK SIZE
MSF-PV2213	100 pcs/pack
MSF-PV4513	100 pcs/pack
MSF-PV2225	100 pcs/pack
MSF-PV4525	100 pcs/pack

**Looking for a different option? We can do it!**  
**Contact us now for more details.**

# PACK SIZES

## Dispenser Dropping Bottles

CATALOG NO.	PACK SIZE
MLS-DDB015	72 pcs/pack
MLS-DDB030	12 pcs/pack
MLS-DDB060	12 pcs/pack
MLS-DDB125	12 pcs/pack
MLS-DDB250	12 pcs/pack

## Aspirator Bottles

CATALOG NO.	PACK SIZE
MLS-34101	1 pc/pack
MLS-34102	1 pc/pack
MLS-34103	1 pc/pack

## Drosophila Bottles

CATALOG NO.	PACK SIZE
MLS-DB-6OZ	200 Pcs/Pack

## Spray Bottles

CATALOG NO.	PACK SIZE
MLS-SB-8OZ	12 pcs/pack
MLS-SB-16OZ	8 pcs/pack

## Wash Bottles

CATALOG NO.	PACK SIZE
MLS-370062	6 pcs/pack
MLS-370061	6 pcs/pack
MLS-370060	6 pcs/pack
MLS-370059	6 pcs/pack
MLS-370058	6 pcs/pack
MLS-370057	6 pcs/pack
MLS-370056	6 pcs/pack
MLS-370055	6 pcs/pack
MLS-370054	6 pcs/pack
MLS-370053	6 pcs/pack
MLS-370052	6 pcs/pack
MLS-370051	6 pcs/pack

## Dropper Bottle

CATALOG NO.	PACK SIZE
MLS-PDB015	24 pcs/pack
MLS-22101	12 pcs/pack
MLS-22102	12 pcs/pack

## Centrifuge Bottles

CATALOG NO.	PACK SIZE
MLS-6330	36 pcs/pack
MLS-6332	24 pcs/pack

## Reagent Bottles

CATALOG NO.	PACK SIZE
MLS-33461	72 pcs/pack
MLS-33462	12 pcs/pack
MLS-33463	12 pcs/pack
MLS-33464	12 pcs/pack
MLS-33465	12 pcs/pack
MLS-33466	6 pcs/pack
MLS-33421	72 pcs/pack
MLS-33422	72 pcs/pack
MLS-33423	72 pcs/pack
MLS-33424	72 pcs/pack
MLS-33425	72 pcs/pack
MLS-33426	12 pcs/pack
MLS-33427	12 pcs/pack
MLS-33428	12 pcs/pack
MLS-33429	6 pcs/pack
MLS-BNMP2L	1 pc/pack
MLS-BNMP4L	1 pc/pack
MLS-BNMP8L	1 pc/pack
MLS-BNM090	500 pcs/pack

Looking for a different option? We can do it!  
Contact us now for more details.

# PACK SIZES

## Media/Storage Bottles Cylindrical Stir Bars

CATALOG NO.	PACK SIZE
MLS-BMPE250	12 pcs/pack
MLS-BMPE500	12 pcs/pack
MLS-BMPE1000	12 pcs/pack
MLS-BMPC150	24 pcs/pack
MLS-BMPC250	24 pcs/pack
MLS-BMPC500	24 pcs/pack
MLS-BMPC1000	24 pcs/pack

## Pivot Ring Stir Bars

CATALOG NO.	PACK SIZE
MLS-SWN501	2 pcs/pack
MLS-SWN502	2 pcs/pack
MLS-SWN503	2 pcs/pack
MLS-SWN504	2 pcs/pack
MLS-SWN505	2 pcs/pack
MLS-SWN506	2 pcs/pack
MLS-SWN507	2 pcs/pack
MLS-SWN508	2 pcs/pack
MLS-SWN509	2 pcs/pack
MLS-SWN510	2 pcs/pack

## Cylindrical Stir Bars

CATALOG NO.	PACK SIZE
MLS-SWN652	2 pcs/pack
MLS-SWN654	2 pcs/pack
MLS-SWN656	2 pcs/pack
MLS-SWN658	2 pcs/pack
MLS-SWN660	2 pcs/pack
MLS-SWN662	2 pcs/pack
MLS-SWN664	2 pcs/pack
MLS-SWN666	2 pcs/pack
MLS-SWN668	2 pcs/pack
MLS-SWN670	2 pcs/pack
MLS-SWN672	2 pcs/pack
MLS-SWN674	2 pcs/pack

## Micro Stir Bars

CATALOG NO.	PACK SIZE
MLS-SWN950	2 pcs/pack
MLS-SWN954	2 pcs/pack
MLS-SWN956	2 pcs/pack

## Oval Stir Bars

CATALOG NO.	PACK SIZE
MLS-SWN520	2 pcs/pack
MLS-SWN522	2 pcs/pack
MLS-SWN524	2 pcs/pack
MLS-SWN525	2 pcs/pack
MLS-SWN526	2 pcs/pack
MLS-SWN528	2 pcs/pack

## Hexagonal Weighing Boats

CATALOG NO.	PACK SIZE
MLS-3615	16 pcs/pack
MLS-3616	16 pcs/pack
MLS-3617	8 pcs/pack
MLS-3618	8 pcs/pack

## Diamond Shaped Weighing Boats

CATALOG NO.	PACK SIZE
MLS-3614-100	12 pcs/pack
MLS-3614-30	16 pcs/pack
MLS-3614-5	36 pcs/pack

**Looking for a different option? We can do it!**  
**Contact us now for more details.**

# PACK SIZES

## Square Weighing Boats

CATALOG NO.	PACK SIZE
MLS-3619-10	24 pcs/pack
MLS-3619-100	12 pcs/pack
MLS-3619-250	6 pcs/pack
MLS-3619-10BK	24 pcs/pack
MLS-3619-100BK	12 pcs/pack
MLS-3619-250BK	6 pcs/pack
MLS-3620	24 pcs/pack
MLS-3621	12 pcs/pack
MLS-3622	6 pcs/pack

## Round Aluminum Dishes with Tabs

CATALOG NO.	PACK SIZE
MLS-D28-500	MLS-D28-500
MLS-D43-100	MLS-D43-100
MLS-D57-100	MLS-D57-100
MLS-D70-100	MLS-D70-100

## Aluminum Weighing or Drying Pans

CATALOG NO.	PACK SIZE
MLS-D123-50	50 pcs/pack
MLS-D126-50	50 pcs/pack
MLS-D125-100	100 pcs/pack
MLS-D11-25	25 pcs/pack
MLS-D12-25	25 pcs/pack

## Micro Aluminum Weighing Dishes

CATALOG NO.	PACK SIZE
MLS-MAWD-0400	250 pcs/pack
MLS-MAWD-1100	250 pcs/pack
MLS-MAWD-1200	250 pcs/pack
MLS-MAWD-1201	250 pcs/pack
MLS-MAWD-1120	50 pcs/pack
MLS-MAWD-3500	100 pcs/pack
MLS-MAWD-0350	250 pcs/pack

## General Purpose Aluminum Weighing Dishes

CATALOG NO.	PACK SIZE
MLS-D44-100	100 pcs/pack
MLS-D63-100	100 pcs/pack
MLS-D75-100	100 pcs/pack
MLS-D200-50	50 pcs/pack
MLS-D500-50	50 pcs/pack
MLS-111	1000 pcs/pack

## Smooth-Walled Aluminum Weighing Dishes

CATALOG NO.	PACK SIZE
MLS-D57S-100	100 pcs/pack
MLS-D70S-100	100 pcs/pack
MLS-D140S-100	100 pcs/pack

**Looking for a different option? We can do it!**  
**Contact us now for more details.**