



Laboratory Consumables and Supplies

HEALTH AND SAFETY PRODUCTS



CATALOG NO.	DESCRIPTION
MLS-GOGGLE-P	Safety Goggles
Flexible safety goggles are manufactured from PVC and are designed to protect against splashes and impact. Feature four vents to prevent fogging and an adjustable band for a comfortable fit. Safety goggles are ANSI Z87.1-2020 certified.	

CATALOG NO.	DESCRIPTION
MLS-HNDPTC	Hand Protector, Silicone Rubber
Made of silicone rubber and features special grips to prevent containers from slipping. Ideal for gripping hot objects such as beakers, flasks and bottles.	



CATALOG NO.	DESCRIPTION
MLS-P70106	Safety Goggles Box, Acrylic, Small
MLS-P70107	Safety Goggles Box, Acrylic, Large
Store goggles, glasses, or other small lab products using this clear acrylic box. It can be used on a desktop or can be mounted on a wall. Includes a lid to keep the box free of dust or contamination.	

CATALOG NO.	DESCRIPTION
MLS-SPB-01	Plastic bandages. Promote healing and offer protection against infection



CATALOG NO.	DESCRIPTION
MLS-SAB-02	Adhesive Bandage Dispenser Refill

HEALTH AND SAFETY PRODUCTS



CATALOG NO.	DESCRIPTION
MLS-NGXS-1000	Examination Grade Gloves, X-Small, Nitrile, 4.5-mil, Powder-Free, Blue
MLS-NGS-1000	Examination Grade Gloves, Small, Nitrile, 4.5-mil, Powder-Free, Blue
MLS-NGM-1000	Examination Grade Gloves, Medium, Nitrile, 4.5-mil, Powder-Free, Blue
MLS-NGL-1000	Examination Grade Gloves, Large, Nitrile, 4.5-mil, Powder-Free, Blue
MLS-NGXL-1000	Examination Grade Gloves, X-Large, Nitrile, 4.5-mil, Powder-Free, Blue
MLS-NG2XL-1000	Examination Grade Gloves, 2X-Large, Nitrile, 4.5-mil, Powder-Free, Blue
Pure 100% premium nitrile; Textured surface provides superior wet and dry grip; Exceptional tactile sensitivity; Beaded cuff for added strength; Protects against oil, grease and organic solvents; Superior puncture resistance; Latex-free	



CATALOG NO.	DESCRIPTION
MLS-SK-03	Kimtech™ Science Kimwipes™ Delicate Task Wipes, Specialty. Available in different sizes.

CATALOG NO.	DESCRIPTION
MLS-SFM-04	Disposable Face Masks with Adjustable nose clip & ear loops provide a custom fit.



CATALOG NO.	DESCRIPTION
MLS-SB-05	First Aid Blanket, Multi-Blend Fibre, Medical Device Class 1

HEALTH AND SAFETY PRODUCTS



CATALOG NO.	DESCRIPTION
MLS-LCMXS1	Laboratory Coats, Poly- Cotton, Extra Small, Men, Size 34
MLS-LCMS01	Laboratory Coats, Poly- Cotton, Small, Men, Size 36
MLS-LCMM01	Laboratory Coats, Poly- Cotton, Medium, Men, Size 40
MLS-LCML01	Laboratory Coats, Poly- Cotton, Large, Men, Size 44
MLS-LCMXL1	Laboratory Coats, Poly- Cotton, Extra Large, Men, Size 48
MLS-LCMXXL1	Laboratory Coats, Poly- Cotton, Double Extra Large, Men, Size 50
MLS-LCM3XL1	Laboratory Coats, Poly- Cotton, Triple Extra Large, Men, Size 52
MLS-LCLS01	Laboratory Coats, Poly- Cotton, Small, women, Size 4-6
MLS-LCLM01	Laboratory Coats, Poly- Cotton, Medium, women, Size 8-10
MLS-LCLL01	Laboratory Coats, Poly- Cotton, Large, women, Size 12-14
MLS-LCLXL1	Laboratory Coats, Poly- Cotton, Extra large, women, Size 16-18
MLS-LCLXXL1	Laboratory Coats, Poly- Cotton, Double Extra Large, women, Size 20-22

CATALOG NO.	DESCRIPTION
MLS-SPP-S	Lab Coat, Polypropylene, White, Small
MLS-SPP-M	Lab Coat, Polypropylene, White, Medium
MLS-SPP-L	Lab Coat, Polypropylene, White, Large
MLS-SPP-XL	Lab Coat, Polypropylene, White, X-Large
MLS-SPP-2X	Lab Coat, Polypropylene, White, 2X-Large
MLS-SPP-3X	Lab Coat, Polypropylene, White, 3X-Large
MLS-SPP-4X	Lab Coat, Polypropylene, White, 4X-Large



HEALTH AND SAFETY PRODUCTS



CATALOG NO.	DESCRIPTION
MLS-SPC-40	Lab Coats, Poly-Cotton, 40, Navy Blue
MLS-SPC-42	Lab Coats, Poly-Cotton, 42, Navy Blue
MLS-SPC-44	Lab Coats, Poly-Cotton, 44, Navy Blue

CATALOG NO.	DESCRIPTION
MLS-LBCTSM	Tie-dyed Laboratory Coats, Poly- Cotton, Small
MLS-LBCTMD	Tie-dyed Laboratory Coats, Poly- Cotton, Medium
MLS-LBCTLG	Tie-dyed Laboratory Coats, Poly- Cotton, Large
MLS-LBCTXL	Tie-dyed Laboratory Coats, Poly- Cotton, Extra Large
MLS-LBCT2XL	Tie-dyed Laboratory Coats, Poly- Cotton, Double Extra large



CATALOG NO.	DESCRIPTION
MLS-CNS428-S	Protective Coverall with Hood, ProShield® 50 microporous film, Small
MLS-CNS428-M	Protective Coverall with Hood, ProShield® 50 microporous film, Medium
MLS-CNS428-L	Protective Coverall with Hood, ProShield® 50 microporous film, Large
MLS-CNS428-XL	Protective Coverall with Hood, ProShield® 50 microporous film, X-Large
MLS-CNS428-2XL	Protective Coverall with Hood, ProShield® 50 microporous film, 2X-Large
MLS-CNS428-3XL	Protective Coverall with Hood, ProShield® 50 microporous film, 3X-Large
MLS-CNS428-4XL	Protective Coverall with Hood, ProShield® 50 microporous film, 4X-Large
MLS-CNS428-5XL	Protective Coverall with Hood, ProShield® 50 microporous film, 5X-Large
MLS-CNS428-6XL	Protective Coverall with Hood, ProShield® 50 microporous film, 6X-Large

HEALTH AND SAFETY PRODUCTS



CATALOG NO.	DESCRIPTION
MLS-SFA-KIT	First Aid Kit, Class 1 Medical Device, Plastic Box

CATALOG NO.	DESCRIPTION
MLS-SSFA-KIT	Office Standard First Aid Kits, Class 1 Medical Device, Plastic Box



CATALOG NO.	DESCRIPTION
MLS-SPC-05	Pleated Bouffant Cap, Polypropylene, Blue. Available in different sizes.

CATALOG NO.	DESCRIPTION
MLS-SBN-06	Beard Nets, Polypropylene, White



CATALOG NO.	DESCRIPTION
MLS-SSC-07	Shoe Covers, Polypropylene, Blue. Available in different sizes

CATALOG NO.	DESCRIPTION
MLS-SSP-08	KleenGuard™ Sleeve Protector, 18" long, Microporous, White



HEALTH AND SAFETY PRODUCTS



CATALOG NO.	DESCRIPTION
MLS-SDT-09	Jumbo Thermometer with Fever Glow, Digital

CATALOG NO.	DESCRIPTION
MLS-SHFP-10	Head and Face Protection Systems



CATALOG NO.	DESCRIPTION
MLS-SIC-11	Instant Compress, Cold, Single Use, 5" x 6"

CATALOG NO.	DESCRIPTION
MLS-SIC-12	Instant Compress, Hot, Single Use, 6" x 9"



CATALOG NO.	DESCRIPTION
MLS-SHC-KIT	Hand Cleansing Kit, Towelette

CATALOG NO.	DESCRIPTION
MLS-SHW-13	Antimicrobial Hand Wipes, Towelette, Antibiotic



HEALTH AND SAFETY PRODUCTS



CATALOG NO.	DESCRIPTION
MLS-SSD-14	Quatromycide V Disinfectant, Sanitizer and Deodorizer, Jug



CATALOG NO.	DESCRIPTION
MLS-SSC-15	Omni-Cleaner™ XL is a steam autoclave cleaner.



CATALOG NO.	DESCRIPTION
MLS-SLA-16	Topical First Aid Treatment, Liquid, Antiseptic



CATALOG NO.	DESCRIPTION
MLS-STT-17	Proviodyne Topical Treatment, Liquid, Antiseptic



CATALOG NO.	DESCRIPTION
MLS-SBO-18	Bacitin Topical Treatment, Ointment



CATALOG NO.	DESCRIPTION
MLS-SSL-19	Instrument Sterilizing Solution, Liquid

LABORATORY CONSUMABLES

Membra-cel® Dialysis membranes

Seamless cellulose membranes featuring sharp Molecular Weight Cut-Offs in both tubing and flat sheets. MEMBRA-CEL® dialysis membranes are made of regenerated cellulose using the viscose process. These semi-permeable membranes have distinct Molecular Weight Cut-Off (MWCOs) of 14,000, 7,000 and 3,500 daltons by excluding macromolecules with molecular radii from 50 down to 35 angstroms.



CATALOG NO.	MOLECULAR WEIGHT CUT-OFF (KILODALTONS)	FLAT WIDTH	DIAMETER	THICKNESS
MLS-MC-MD10	3.5, 7, and 14	10mm	6mm	51µm
MLS-MC-MD25	3.5, 7, and 14	25mm	16mm	20µm
MLS-MC-MD34	3.5, 7, and 14	34mm	22mm	23µm
MLS-MC-MD44	3.5, 7, and 14	44mm	28mm	20µm
MLS-MC-MD55	3.5, 7, and 14	55mm	35mm	25µm
MLS-MC-MD77	3.5, 7, and 14	77mm	49mm	41µm

Features include:

- Excellent pH and chemical stability, the membranes are stable in weak acids and bases
- Great compatibility with many salts –CaCl₂, (NH₄)₂SO₄, etc.
- Excellent compatibility with most aqueous organic solvents, such as alcohols and ketones
- Excellent temperature resistance Sensitive to microbial activity when humidified, however, growth of cellulolytic micro-organisms can be prevented by keeping the membrane in water containing benzoate, benzoic acid, or other preservatives
- Very low protein absorption
- Available in widths ranging from 10mm to 172mm and in rolls lengths of 100 or 500 feet
- Contains low levels of contaminants, some membranes contain glycerin as a humectant, contaminants may be removed by boiling prior to use

LABORATORY CONSUMABLES



CATALOG NO.	MLS- 53223-7601
DESCRIPTION	Porcelain Buchner Funnels
CAPACITY	160 ml
PLATE DIAM	130 mm
FILTER PAPER SIZE	12.5 cm
They possess the properties such as resistant to Thermal shocks and Chemical actions.	

CATALOG NO.	DESCRIPTION
MLS-53121-1110	70mL Mortar
MLS-53121-1220	160mL Mortar
MLS-53121-1320	400mL Mortar
MLS-53121-1500	650mL Mortar
MLS-53121-3400	175mm Pestle



CATALOG NO.	DESCRIPTION
MLS-CFT-01	Dial Thermometer Clips. It is used to attach thermometers to beakers, containers, tanks, etc.

CATALOG NO.	DESCRIPTION	THERMOMETER LENGTH
83025A	Thermometer Cases	305 mm
83025E	Thermometer Cases	229 mm
83025H	Thermometer Cases	406 mm
83025N	Thermometer Cases	393 mm



LABORATORY CONSUMABLES



CATALOG NO.	DESCRIPTION	LENGTH
MLS-66CS-8	PVC hydrometer cases	8.5 in.
MLS-66CS-8F	PVC hydrometer cases	8.5 in.
MLS-66CS-13	PVC hydrometer cases	13.5 in.
MLS-66CS-13J	PVC hydrometer cases	13.5 in.
MLS-66CS-16	PVC hydrometer cases	16.5 in.

CATALOG NO.	DESCRIPTION	TYPE	BRIX VALUE (20°C)	R.I. VALUE (20°C)
MLS-45000-OS025	Refractometer, Standards	Oil	60.00%	1.4420 nD
MLS-45000-OS200	Refractometer, Standards	Water	0.00%	1.33297 nD



CATALOG NO.	DESCRIPTION
MLS-LPR100	Red Litmus Paper, Vial 100/Strips
MLS-LPB100	Blue Litmus Paper, Vial 100/Strips
MLS-LPN100	Neutral Litmus Paper, Vial 100/Strips

CATALOG NO.	DESCRIPTION
MLS-FSLT02	Plastic 6" long flashlight includes long-lasting LED illumination and lanyard.



LABORATORY CONSUMABLES



CATALOG NO.	DESCRIPTION
MLS-62100	Microscope Slide Mailer, PP, For 1 slide
MLS-62102	Microscope Slide Mailer, PP, For two slides
MLS-62103	Microscope Slide Mailer, PP, For 5 slides
MLS-62111	Microscope Slide Mailer, PP, For 11 slides

CATALOG NO.	DESCRIPTION	SIZE
MLS-FPR007	Filter Papers, Circular, Grade 1, Medium porosity and flow rate	7 cm
MLS-FPR009	Filter Papers, Circular, Grade 1, Medium porosity and flow rate	9 cm
MLS-FPR011	Filter Papers, Circular, Grade 1, Medium porosity and flow rate	11 cm
MLS-FPR125	Filter Papers, Circular, Grade 1, Medium porosity and flow rate	12.5 cm
MLS-FPR150	Filter Papers, Circular, Grade 1, Medium porosity and flow rate	15 cm
MLS-FPR180	Filter Papers, Circular, Grade 1, Medium porosity and flow rate	18 cm
MLS-FPR240	Filter Papers, Circular, Grade 1, Medium porosity and flow rate	24 cm
MLS-FPR320	Filter Papers, Circular, Grade 1, Medium porosity and flow rate	32 cm



CATALOG NO.	DESCRIPTION
MLS-FLTR01	Rubber Filter Aid. Inside diameter of the socket at top is 32mm, and outside diameter of the disc is 76mm. The bottom of the disc is flat so that it seals tightly against the lip of the flask.

LABORATORY CONSUMABLES

Clear 2 mil polypropylene Autoclave Bags (No Print)

Intended for use in autoclaving non hazardous materials and waste. Polypropylene bags are recommended for maximum temperature exposure of 135° C (275° F) during cycle. Twist ties included allow proper closure for bag expansion and steam venting during sterilization. Bags meet the strength requirements of the 165 gram dart drop standard as performed per ASTM D1709 test. Twist ties included for proper closure.



CATALOG NO.	DIMENSIONS
MLS-1102-0812	8 in x 12 in
MLS-1102-1015	10 in x 15 in
MLS-1102-1224	12 in x 24 in
MLS-1102-2436	24 in x 36 in

CATALOG NO.	DIMENSIONS
MLS-1102-1924	19 in x 24 in
MLS-1102-2436	24 in x 36 in
MLS-1102-3138	31 in x 38 in
MLS-1102-3748	37 in x 48 in

Clear Autoclavable Biohazard Bags

The ideal bag for processing biohazardous waste in an autoclave prior to disposal. Bags meet the strength requirements of 165 gram dart drop standard as performed per ASTM D1709 test. Polypropylene bags are recommended for maximum temperature exposure of 135° C (275° F) and have a temperature indicator patch printed next to Biohazard symbol. Twist ties included for proper closure.

CATALOG NO.	DIMENSIONS
MLS-1112-0812	8 in x 12 in
MLS-1112-1224	12 in x 24 in
MLS-1112-2436	24 in x 36 in
MLS-1112-1924	19 in x 24 in

CATALOG NO.	DIMENSIONS
MLS-1112-2436	24 in x 36 in
MLS-1112-3138	31 in x 38 in
MLS-1112-3748	37 in x 48 in



LABORATORY CONSUMABLES

Red Autoclavable Biohazard Bags

Engineered to meet the most demanding regulatory specifications worldwide. These bags are specifically designed for use in Biosafety Level 3 & 4 (BSL 3 & 4) research labs for handling of very dangerous materials. These Biohazard bags meet all State & Federal requirements for hazardous waste handling. Including Federal DOT Regulation 49 CFR 173.197 for HazMat handling and markings. They exceed the ASTM D1922 Tear Resistance Test of 480 grams in both directions, As well as ASTM 1709 Impact Resistance of 165 grams



CATALOG NO.	DIMENSIONS
MLS-1212-0812	8 in x 12 in
MLS-1212-1224	12 in x 24 in
MLS-1212-1419	14 in x 19 in
MLS-1212-1924	19 in x 24 in
MLS-1212-2430	24 in x 30 in

CATALOG NO.	DIMENSIONS
MLS-1212-2436	24 in x 36 in
MLS-1212-2535	25 in x 35 in
MLS-1212-3138	31 in x 38 in
MLS-1212-3748	37 in x 48 in

Orange Autoclavable Biohazard Bags

The ideal bag for processing Biohazardous waste in an autoclave prior to disposal. Bags meet the strength requirements of 165 gram dart drop standard as performed per ASTM D1709 test. Polypropylene bags are recommended for maximum temperature exposure of 135° C (275° F) and have a temperature indicator patch printed next to Biohazard symbol. Twist ties included for proper closure.

CATALOG NO.	DIMENSIONS
MLS-1312-0812	8 in x 12 in
MLS-1312-1224	12 in x 24 in
MLS-1312-1924	19 in x 24 in

CATALOG NO.	DIMENSIONS
MLS-1312-2436	24 in x 36 in
MLS-1312-2535	25 in x 35 in
MLS-1312-3748	37 in x 48 in



LABORATORY CONSUMABLES

Yellow Autoclavable Biohazard Bags

The ideal bag for processing Biohazardous waste in an autoclave prior to disposal. Bags meet the strength requirements of 165 gram dart drop standard as performed per ASTM D1709 test. Polypropylene bags are recommended for maximum temperature exposure of 135° C (275° F) and have a temperature indicator patch printed next to Biohazard symbol. Twist ties included for proper closure.



CATALOG NO.	DIMENSIONS
MLS-1412-0812	8 in x 12 in
MLS-1412-1419	14 in x 19 in
MLS-1412-1924	19 in x 24 in
MLS-1412-2436	24 in x 36 in
MLS-1412-2836	28 in x 36 in
MLS-1412-3138	21 in x 38 in
MLS-1412-3748	37 in x 48 in

Infecon Bags

The Infecon® specimen transport bag is the perfect choice for the secure transport and handling of diagnostic specimens. Infecon bags are specially manufactured using a patented process which guarantees the elimination of leaks. The handy document pouch prevents misplacing vital paperwork.



CATALOG NO.	DESCRIPTION
MLS-4760	Infecon Specimen Bags, Zipper Closure, 6 x 6" with Document Pouch
MLS-4764	Infecon Specimen Bags, Zipper Closure, 12 x 12" with Document Pouch
MLS-4765	Infecon Specimen Bags, Zipper Closure, 16 x 16" with Document Pouch
MLS-4767	Infecon Specimen Bags, Zipper Closure, 24 x 24" with Document Pouch