

A collection of laboratory plastic syringes is shown in a blue-tinted environment. The syringes are arranged in a cluster, with some in sharp focus and others blurred in the background. The syringes have various markings on their barrels, including numbers like 1, 2, 3, 4, and 5, and the unit 'mL'. The overall scene is lit with a cool, blue light, creating a clinical and scientific atmosphere.

# Laboratory Plastic Syringes

# PLASTIC SYRINGES

## Air-Tite Product Syringes

Air-Tite laboratory and industrial syringes are the crossroads where quality meets economy. A basic 3-part syringe with consistent quality is ideal for application which simply require low cost transfer vessel. Made to our exciting specifications, these syringes are latex-free, non-pyrogenic, and toxic free. The design includes quality print and siliconized gaskets for a consistently smooth plunger motion and positive plunger stop. For Veterinary and Laboratory use only. All sizes have a centric tip except for 30ml, which has an eccentric tip. Packaged in individual sterile poly bags.



## Featured Products

- Air-tite, Luer Slip Syringes For Laboratory And Veterinary Use, Sterile
- Air-tite Luer Lock Syringes For Laboratory And Veterinary Use, Sterile
- Air-tite Bulk Syringes For Laboratory And Veterinary Use, Unsterile

## Features Include:

- High quality
- Economical
- Laboratory and Veterinary Use
- Manufactured to Our Specifications
- Latex Free
- Non-toxic
- Pyrogen Free
- Positive Plunger Stop
- Standard 3-Part Syringes
- Lubricated Rubber Tip Plunger

**Looking for a different option? We can do it!**  
**Contact us now for more details.**

# AIR-TITE SYRINGES

## Air-Tite Luer Slip Syringes



CATALOG NO.	DESCRIPTION	PACK SIZE
MSY-0001LSB	1cc Tuberculin	100 pcs/pack
MSY-0003LSB	3cc Centric	100 pcs/pack
MSY-0005LSB	5cc Centric	100 pcs/pack
MSY-00010LSB	10cc Centric	100 pcs/pack
MSY-00020LSB	20cc Eccentric	50 pcs/pack
MSY-00030LSB	30cc Eccentric	50 pcs/pack
MSY-00060LSB	60cc Eccentric	25 pcs/pack

## Air-Tite Luer Lock Syringes

CATALOG NO.	DESCRIPTION	PACK SIZE
MSY-0001LLB	1cc	100 pcs/pack
MSY-0003LLB	3cc	100 pcs/pack
MSY-0005LLB	5cc	100 pcs/pack
MSY-00010LLB	10cc	100 pcs/pack
MSY-00020LLB	20cc	50 pcs/pack
MSY-00030LLB	30cc	50 pcs/pack
MSY-00060LLB	60cc	25 pcs/pack



# AIR-TITE SYRINGES

## Air-Tite Bulk Unsterile Syringes



CATALOG NO.	SIZE	TIP TYPE	PACK SIZE	CASE SIZE
MSY-0001LSBB	1mL	Luer Slip	450 pcs/pack	4500 pcs/case
MSY-0001LLBB	1mL	Luer Lock	450 pcs/pack	4500 pcs/case
MSY-0003LSBB	3mL	Luer Slip	350 pcs/pack	3500 pcs/case
MSY-0003LLBB	3mL	Luer Lock	350 pcs/pack	3500 pcs/case
MSY-0005LSBB	5mL	Luer Slip	200 pcs/pack	2000 pcs/case
MSY-0005LLBB	5mL	Luer Lock	200 pcs/pack	2000 pcs/case
MSY-00010LSBB	10mL	Luer Slip	130 pcs/pack	1300 pcs/case
MSY-00010LLBB	10mL	Luer Lock	130 pcs/pack	1300 pcs/case
MSY-00020LSBB	20mL	Luer Slip	50 pcs/pack	500 pcs/case
MSY-00020LLBB	20mL	Luer Lock	50 pcs/pack	500 pcs/case
MSY-00030LSBB	30mL	Luer Slip	50 pcs/pack	500 pcs/case
MSY-00030LLBB	30mL	Luer Lock	50 pcs/pack	500 pcs/case
MSY-00060LLBB	60mL	Luer Lock	25 pcs/pack	250 pcs/case

# PLASTIC SYRINGES

## Norm-Ject and Hanke-Ject

NORM-JECT<sup>®</sup> and Henke-Ject<sup>®</sup> syringes are latex free, and contain no rubber, silicone oil, styrene or DEHP. These syringes are the choice for any situation needing an inert, non-reactive syringe. Manufactured only from laboratory grade polypropylene and polyethylene. These unique plastic syringes are more chemically resistant than rubber-tipped syringes and have a positive safety stop to help prevent accidental spills. NORM-JECT<sup>®</sup> and Henke-Ject<sup>®</sup> syringes are individually sterile strip packed.

### Featured Products

- B.braun Norm-ject (1ml To 20ml) And Hsw Henke Ject (30ml To 50ml), All-plastic, 2-piece, Luer Slip Syringes, Sterile
- B.braun Norm-ject (3ml To 20ml) And Hsw Henke Ject (30ml To 50ml), All-plastic, 2-piece, Luer Lock Syringes, Sterile
- Hsw Henke-ject, All-plastic, 2-piece, Catheter Tip Syringes, Sterile
- Hsw All-plastic, 2-piece, Laboratory Convenience Packs, Unsterile
- Hsw All-plastic, 2-piece, Bulk Non-sterile

## Hanke-Sass Wolf 3-Part Syringes

HSW HENKE-JECT<sup>®</sup> (Previously known as Soft-Ject<sup>®</sup>) syringes are a high quality disposable rubber piston syringe offering substantial savings. A lubricated synthetic rubber tipped piston allows for smooth draw while a positive stop prevents accidental pullout. These syringes are ETO gas sterilized and are latex and pyrogen free. All HENKE-JECT<sup>®</sup> syringes are packaged in individual sterile strip packs. FDA registered CE Mark, ISO 9001 certified.

### Featured Products

- Hsw Soft-ject / Henke-ject, 3-part, Luer Slip Syringes, Sterile
- Hsw Soft-ject / Henke-ject, 3-part, Luer Lock Syringes, Sterile

**Looking for a different option? We can do it!**  
**Contact us now for more details.**

# NORM-JECT SYRINGES

## Norm-Ject Luer Slip 2-part Syringes



CATALOG NO.	DESCRIPTION	PACK SIZE
MSY-HSWNJLS1	1mL Tuberculin	100 pcs/pack
MSY-HSWNJLS3	3mL Centric	100 pcs/pack
MSY-HSWNJLS5	5mL Centric	100 pcs/pack
MSY-HSWNJLS10	10mL Eccentric	100 pcs/pack
MSY-HSWNJLS20	20mL Eccentric	50 pcs/pack

## Norm-Ject Luer Lock 2-part Syringes

CATALOG NO.	DESCRIPTION	PACK SIZE
MSY-HSWNJLL3	3mL	100 pcs/pack
MSY-HSWNJLL5	5mL	100 pcs/pack
MSY-HSWNJLL10	10mL	100 pcs/pack
MSY-HSWNJLL20	20mL	100 pcs/pack



# HSW HENKE-JECT SYRINGES

## HSW Henke-Ject Luer Slip 2-part Syringes



CATALOG NO.	DESCRIPTION	PACK SIZE
MSY-HSWHJLS30	30mL Eccentric	50 pcs/pack
MSY-HSWHJLS50	50mL Eccentric	30 pcs/pack

# HSW HENKE-JECT SYRINGES

## HSW Henke-Ject Luer Lock 2-part Syringes



CATALOG NO.	DESCRIPTION	PACK SIZE
MSY-HSWHJLL30	30mL	50 pcs/pack
MSY-HSWHJLL50	50mL	30 pcs/pack

## HSW Henke-Ject Catheter Tip Syringes

CATALOG NO.	DESCRIPTION	PACK SIZE
MSY-HSWCT50	50mL	30 pcs/pack



# 2-PART UNSTERILE SYRINGES



CATALOG NO.	SIZE	TIP TYPE	CASE SIZE
MSY-0001B	1mL	Luer Slip	7,000 pcs/case
MSY-0003CB	3mL Centric	Luer Slip	6,300 pcs/case
MSY-0003B	3mL	Luer Lock	6,300 pcs/case
MSY-0005CB	5mL Centric	Luer Slip	3,600 pcs/case
MSY-0005B	5mL	Luer Lock	3,600 pcs/case
MSY-0010EB	10mL Eccentric	Luer Slip	2,000 pcs/case
MSY-0010B	10mL	Luer Lock	2,000 pcs/case
MSY-0020EB	20mL Eccentric	Luer Slip	1,000 pcs/case
MSY-0020B	20mL	Luer Lock	1,000 pcs/case

# SYRINGES WITH NEEDLES

## Air-Tite Luer Slip Syringes With Needles



CATALOG NO.	SIZE	PACK SIZE
MSY-0125LSN	1cc x 25g x 5/8"	100 pcs/pack
MSY-0127LSN	1cc x 27g x 1/2"	100 pcs/pack
MSY-0318LSN	3cc x 18g x 1"	100 pcs/pack
MSY-0320LSN	3cc x 20g x 1"	100 pcs/pack
MSY-0322LSN	3cc x 22g x 3/4"	100 pcs/pack
MSY-0322ILSN	3cc x 22g x 1"	100 pcs/pack
MSY-0325LSN	3cc x 25g x 5/8"	100 pcs/pack

## Air-Tite Luer Lock Syringes with Needles

CATALOG NO.	DESCRIPTION	PACK SIZE
MSY-0125ILLN	1cc x 25g x 5/8"	100 pcs/pack
MSY-0318LLN	3cc x 18g x 1"	100 pcs/pack
MSY-0320LLN	3cc x 20g x 1"	100 pcs/pack
MSY-0320ILLN	3cc x 20g x 1-1/2"	100 pcs/pack
MSY-0322LLN	3cc x 22g x 3/4"	100 pcs/pack
MSY-0322ILLN	3cc x 22g x 1"	100 pcs/pack
MSY-0325LLN	3cc x 25g x 5/8"	100 pcs/pack
MSY-0518LLN	5cc x 18g x 1"	100 pcs/pack
MSY-0320LLN	5cc x 20g x 1-1/2"	100 pcs/pack
MSY-1018LLN	10cc x 18g x 1-1/2"	100 pcs/pack
MSY-2018LLN	20cc x 18g x 1-1/2"	50 pcs/pack





# SYRINGES WITH NEEDLES

## Henke-Ject Syringes with Fixed Needle



CATALOG NO.	SIZE	TIP TYPE	PACK SIZE
MLS-RUI3012	1cc x 30g x 1/2" tuberculin	Fixed	100 pcs/pack

## Henke-Ject Syringes with By-pack Needle



CATALOG NO.	SIZE	TIP TYPE	PACK SIZE
MLS-R22234	2.5mL x 22g x 3/4"	Luer Slip	100 pcs/pack
MLS-R2221	2.5mL x 22g x 1"	Luer Slip	100 pcs/pack
MLS-R22558	2.5mL x 25g x 5/8"	Luer Slip	100 pcs/pack
MLS-R2018	20mL x 18g x 1-1/2"	Luer Slip	100 pcs/pack

## Exel Syringes with mounted needles

CATALOG NO.	SIZE	TIP TYPE	PACK SIZE
MLS-ES1251	1cc x 25g x 1"	Luer Slip	100 pcs/pack
MLS-ES3201	3cc x 20g x 1"	Luer Slip	100 pcs/pack
MLS-ES3221	3cc x 22g x 1"	Luer Slip	100 pcs/pack
MLS-EL3201	3cc x 20g x 1"	Luer Lock	100 pcs/pack
MLS-EL5201	5cc x 20g x 1"	Luer Lock	100 pcs/pack
MLS-EL10201	10cc x 20g x 1"	Luer Lock	100 pcs/pack



# SYRINGE NEEDLES

## Air-Tite Brand Hypodermic Needles



CATALOG NO.	SIZE	WALL TYPE	PACK SIZE
MLS-N142	14g x 2"	Regular	100 pcs/pack
MLS-N161	16g x 1"	Regular	100 pcs/pack
MLS-N181	18g x 1"	Regular	100 pcs/pack
MLS-N193D	19g x 3"	Regular	100 pcs/pack
MLS-N201	20g x 1"	Regular	100 pcs/pack
MLS-N211	21g x 1"	Regular	100 pcs/pack
MLS-N221	22g x 1"	Regular	100 pcs/pack
MLS-N231	23g x 1"	Regular	100 pcs/pack
MLS-N251	25g x 1"	Regular	100 pcs/pack
MLS-N2612	26g x 1/2"	Regular	100 pcs/pack
MLS-N2712	27g x 1/2"	Regular	100 pcs/pack

## Air-Tite Brand Metal Hub Needles

CATALOG NO.	SIZE	WALL TYPE	PACK SIZE
MLS-NMH19M	19g x 1-1/2" non-silicone	Regular	12 pcs/pack
MLS-NMH2014	20g x 1/4"	Regular	10 pcs/pack
MLS-NMH2012	20g x 1/2"	Regular	10 pcs/pack



# SYRINGE NEEDLES

## Henke-Ject Hypodermic Needles



CATALOG NO.	SIZE	WALL TYPE	PACK SIZE
MLS-NH15112	15g x 1-1/2"	Thin	100 pcs/pack
MLS-NH172	17g x 2"	Thin	100 pcs/pack
MLS-NH18112	18g x 1-1/2"	Thin	100 pcs/pack
MLS-NH182	18g x 2"	Thin	100 pcs/pack
MLS-NH201	20g x 1"	Thin	100 pcs/pack
MLS-NH212	21g x 2"	Thin	100 pcs/pack
MLS-NH221	22g x 1"	Thin	100 pcs/pack
MLS-NH222	22g x 2"	Thin	100 pcs/pack
MLS-NH231	23g x 1"	Thin	100 pcs/pack
MLS-NH251	25g x 1"	Thin	100 pcs/pack
MLS-NH272	27g x 2"	Regular	100 pcs/pack
MLS-NH3012	30g x 1/2"	Regular	100 pcs/pack

## HSW Eco Metal Hub Needles

CATALOG NO.	SIZE	WALL TYPE	PACK SIZE
MLS-NHM1634	16g x 3/4"	Regular	12 pcs/pack
MLS-NHM1838	18g x 3/8"	Regular	12 pcs/pack
MLS-NHM2038	20g x 3/8"	Regular	12 pcs/pack
MLS-NHM2238	22g x 3/8"	Regular	12 pcs/pack

